

# **SIMPA**

**version 0.4.0**

**CAMI (Computer Assisted Medical Interventions), DKFZ, Heidelberg and Cancer Research  
UK, Cambridge Institute (CRUK CI)**

June 03, 2021




# Contents

<b>Welcome to the SIMPA documentation!</b>	<b>1</b>
<b>README</b>	<b>1</b>
SIMPA Install Instructions	1
Building the documentation	1
External Tools installation instructions	1
mcx (Optical Forward Model)	1
k-Wave (Acoustic Forward Model)	2
Overview	2
Simulating photoacoustic images	2
Path Management	2
How to contribute	2
Performance profiling	3
<b>Developer Guide</b>	<b>3</b>
How to contribute	3
Coding style	3
Documenting your code	4
Adding literature absorption spectra	4
<b>Class references</b>	<b>4</b>
Module: core	5
Volume creation	5
Model-based volume creation	5
Segmentation-based volume creation	6
Optical forward modelling	6
mcx integration	7
Acoustic forward modelling	7
k-Wave integration	8
Image reconstruction	9
Backprojection	10
Delay-Multiply-And-Sum (DMAS)	11
signed Delay-Multiply-And-Sum (sDMAS)	12
Time Reversal	12
Processing Components	13
Noise Models	13
Digital device twins	15
Detection Geometries	15
Illumination Geometries	18
Models of real world devices	19
Module: utils	21
Libraries	38
Module: io_handling	44

Module: log	44
<b>Examples</b>	<b>45</b>
Performing an optical forward simulation	45
Performing a complete forward simulation with acoustic modeling, optical modeling, as well as image reconstruction	48
Defining custom tissue structures and properties	52
Defining a custom digital device twin class	53
Defining custom tissue types	53
Load a segmentation mask and use it to simulate	54
<b>Index</b>	<b>57</b>
<b>Python Module Index</b>	<b>67</b>

## Welcome to the SIMPA documentation!



## README

The Simulation and Image Processing for Photoacoustic Imaging (SIMPA) toolkit.

## SIMPA Install Instructions

You can install simpa with pip. Simply run:

```
pip install simpa
```

For a manual installation from the code, please follow steps 1 - 3:

1. `git clone https://github.com/CAMI-DKFZ/simpa.git`
2. `git checkout master`
3. `git pull`

Now open a python instance in the 'simpa' folder that you have just downloaded. Make sure that you have your preferred virtual environment activated

1. `cd simpa`
2. `pip install -r requirements.txt`
3. `python -m setup.py install`
4. Test if the installation worked by using `python` followed by `import simpa` then `exit()`

If no error messages arise, you are now setup to use simpa in your project.

You also need to manually install the pytorch library to use all features of SIMPA. To this end, use the pytorch website tool to figure out which version to install: <https://pytorch.org/get-started/locally/>

## Building the documentation

When the installation went fine and you want to make sure that you have the latest documentation you should do the following steps in a command line:

1. Navigate to the `simpa` source directory (same level where the `setup.py` is in)
2. Execute `sphinx-build -b pdf -a simpa_documentation/src simpa_documentation` the command
3. Find the PDF file in `simpa_documentation/simpa_documentation.pdf`

## External Tools installation instructions

### *mcx (Optical Forward Model)*

Either download suitable executables or build yourself from the following sources:

<http://mcx.space/>

In order to obtain access to all custom sources that we implemented, please build mcx yourself from the following mcx Github fork: <https://www.github.com/jgroehl/mcx>

For the installation, please follow the instructions from the original repository. Please note that there might be compatibility issues using mcx-cl with the MCX Adapter as this use case is not being tested and supported by the SIMPA developers.

## ***k-Wave (Acoustic Forward Model)***

Please follow the following steps and use the k-Wave install instructions for further (and much better) guidance under:

<http://www.k-wave.org/>

1. Install MATLAB with the core and parallel computing toolboxes activated at the minimum.
2. Download the kWave toolbox
3. Add the kWave toolbox base bath to the toolbox paths in MATLAB
4. If wanted: Download the CPP and CUDA binary files and place them in the k-Wave/binaries folder
5. Note down the system path to the `matlab` executable file.

## **Overview**

The main use case for the simpa framework is the simulation of photoacoustic images. However, it can also be used for image processing.

## ***Simulating photoacoustic images***

A basic example on how to use simpa in your project to run an optical forward simulation is given in the `samples/minimal_optical_simulation.py` file.

## ***Path Management***

As a pipelining tool that serves as a communication layer between different numerical forward models and processing tools, SIMPA needs to be configured with the paths to these tools on your local hard drive. To this end, we have implemented the `PathManager` class that you can import to your project using `from simpa.utils import PathManager`. The `PathManager` looks for a `path_config.env` file (just like the one we provided in the `simpa_examples`) in the following places in this order:

1. The optional path you give the `PathManager`
2. Your `$HOME$` directory
3. The current working directory
4. The SIMPA home directory path

## **How to contribute**

Please find a more detailed description of how to contribute as well as code style references in our **`developer_guide.md`**

The SIMPA code is written and maintained on a closed git repository that is hosted on a server of the German Cancer Research Center (DKFZ), Heidelberg, Germany and changes to the `develop` or `master` branch are mirrored on Github. As such, only the current `master` and `develop` branch of the repository are open source.

To contribute to SIMPA, please fork the SIMPA github repository and create a pull request with a branch containing your suggested changes. The core team developers will then review the suggested changes and integrate these into the code base.

Please make sure that you have included unit tests for your code and that all previous tests still run through.

There is a regular SIMPA status meeting every Friday on even calendar weeks at 10:00 CET/CEST and you are very welcome to participate and raise any issues or suggest new features. You can join the meeting using the following link:

<https://meet.google.com/rze-bxej-cvj>

Please see the github guidelines for creating pull requests:  
<https://docs.github.com/en/github/collaborating-with-issues-and-pull-requests/about-pull-requests>

## Performance profiling

Do you wish to know which parts of the simulation pipeline cost the most amount of time? If that is the case then you can use the following commands to profile the execution of your simulation script. You simply need to replace the `myscript` name with your script name.

```
python -m cProfile -o myscript.cprof myscript.py
```

```
pyprof2calltree -k -i myscript.cprof
```

## Developer Guide

Dear SIMPA developers, Dear person who wants to contribute to the SIMPA toolkit,

First of all: Thank you for your participation and help! It is much appreciated! This Guide is meant to be used as a collection of How-To's to contribute to the framework. In case you have any questions, do not hesitate to get in touch with the members of the core development team:

Kris K. Dreher ([k.dreher@dkfz-heidelberg.de](mailto:k.dreher@dkfz-heidelberg.de))

Jane M. Groehl ([jane.groehl@cruk.cam.ac.uk](mailto:jane.groehl@cruk.cam.ac.uk))

## How to contribute

The SIMPA code is written and maintained on a closed repository that is hosted on a server of the German Cancer Research Center and changes to the `develop` or `master` branch are mirrored on Github (<https://github.com/CAMI-DKFZ/simpa/>). As such, only the current `master` and `develop` branch of the repository are open source.

To make us aware of an issue, please create an issue on the SIMPA github repository.

To contribute to SIMPA, please fork the SIMPA github repository and create a pull request with a branch containing your suggested changes. The core team developers will then review the suggested changes and integrate these into the code base.

Please make sure that you have included unit tests for your code and that all previous tests still run through.

There is a regular SIMPA status meeting every Friday on even calendar weeks at 10:00 CET/CEST and you are very welcome to participate and raise any issues or suggest new features. You can obtain the meeting links from the core developer team. We also have a Slack workspace that you can join if you are interested to contribute.

Please see the github guidelines for creating pull requests:  
<https://docs.github.com/en/github/collaborating-with-issues-and-pull-requests/about-pull-requests>

## Coding style

When writing code for SIMPA, please use the PEP 8 python coding conventions (<https://www.python.org/dev/peps/pep-0008/>) and consider using the following structures in your code in order to make a new developer or someone external always know exactly what to expect.

- Classnames are written in camel-case notation `ClassName`
- Function names are written in small letter with `_` as the delimiter `function_name`
- Function parameters are always annotated with their type `arg1: type = default`

- Only use primitive types as defaults. If a non-primitive type is used, then the default should be `None` and the parameter should be initialized in the beginning of a function.
- A single line of code should not be longer than 120 characters.
- Functions should follow the following simple structure:
  1. Input validation (arguments all not `None`, correct type, and acceptable value ranges?)
  2. Processing (clean handling of errors that might occur)
  3. Output generation (sanity checking of the output before handing it off to the caller)

## Documenting your code

Only documented code will appear in the sphinx generated documentation.

A class should be documented using the following syntax:

```
class ClassName(Superclass):  
    """  
    Explain how the class is used and what it does.  
    """
```

For functions, a lot of extra attributes can be added to the documentation:

```
def function_name(self, arg1:type = default, arg2:type = default) -> return_type:  
    """  
    Explain how the function is used and what it does.  
  
    :param arg1: type, value range, Null acceptable?  
    :param arg2: type, value range, Null acceptable?  
    :returns: type, value range, does it return Null?  
    :raises ExceptionType: explain when and why this exception is raised  
    """
```

## Adding literature absorption spectra

The central point, where absorption spectra are collected and handled is in `simpa.utils.libraries.spectra_library.py`. The file comprises the class `AbsorptionSpectrumLibrary`, in which the new absorption spectra can be added using the following two steps:

1. In the beginning of the class, there is a bunch of constants that define spectra using the `AbsorptionSpectrum` class. Add a new constant here:  
`NEW_SPECTRUM = AbsorptionSpectrum(absorber_name, wavelengths, absorptions)`. By convention, the naming of the constant should be the same as the `absorber_name` field. The `wavelengths` and `absorptions` arrays must be of the same length and contain corresponding values.
2. In the `__init__` method of the `AbsorptionSpectrumLibrary` class, the class constants are added to an internal list. This has the benefit of enabling the Library class to be iterable. Add your newly added constant field to the list here.
3. Your absorption spectrum is now usable throughout all of `simpa` and is accessible using the `SPECTRAL_LIBRARY` singleton that can be imported using `from simpa.utils import SPECTRAL_LIBRARY`.

## Class references

This description details the three principle modules of the SIMPA toolkit and gives an insight into their constituents. The core is concerned with providing interfaces for the simulation tools, while the `utils` module contains many scripts and classes to facilitate the use of the simulation pipeline.



## Module: core

The purpose of the core module is to provide interfaces that facilitate the integration of toolboxes and code for photoacoustic modeling into a single continuous pipeline.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
simpa.core.simulation.simulate(simulation_pipeline: list, settings:
simpa.utils.settings.Settings, digital_device_twin:
simpa.core.device_digital_twins.digital_device_twin_base.DigitalDeviceTwinBase)
```

This method constitutes the starting point for the simulation pipeline of the SIMPA toolkit.

### Parameters:

- **simulation\_pipeline** – a list of callable functions
- **settings** – settings dictionary containing the simulation instructions
- **digital\_device\_twin** – a digital device twin of an imaging device as specified by the DigitalDeviceTwinBase class.

### Raises:

- **TypeError** – if one of the given parameters is not of the correct type
- **AssertionError** – if the digital device twin is not able to simulate the settings specification

**Returns:** list with the save paths of the simulated data within the HDF5 file.

## Volume creation

The core contribution of the SIMPA toolkit is the creation of in silico tissue-mimicking phantoms. This feature is represented by the volume\_creation module, that two main volume creation modules:

- Model-based creation of volumes using a set of rules
- Segmentation-based creation of volumes

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.volume_creation_module.VolumeCreatorModuleBase (global_settings:
simpa.utils.settings.Settings)
```

Use this class to define your own volume creation adapter.

**abstract create\_simulation\_volume ()** → dict

This method will be called to create a simulation volume.

**run (device)**

Executes the respective simulation module

**Parameters:** **digital\_device\_twin** – The digital twin that can be used by the digital device\_twin.

## Model-based volume creation

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.volume_creation_module.volume_creation_module_model_based_adapter.VolumeCreationModelModelBasedAdapter (global_settings: simpa.utils.settings.Settings)
```

The model-based volume creator uses a set of rules how to generate structures to create a simulation volume. These structures are added to the dictionary and later combined by the algorithm:

```
# Initialise settings dictionaries
simulation_settings = Settings()
all_structures = Settings()
structure = Settings()
```

```

# Definition of an example structure.
# The concrete structure parameters will change depending on the
# structure type
structure[Tags.PRIORITY] = 1
structure[Tags.STRUCTURE_START_MM] = [0, 0, 0]
structure[Tags.STRUCTURE_END_MM] = [0, 0, 100]
structure[Tags.MOLECULE_COMPOSITION] = TISSUE_LIBRARY.muscle()
structure[Tags.CONSIDER_PARTIAL_VOLUME] = True
structure[Tags.ADHERE_TO_DEFORMATION] = True
structure[Tags.STRUCTURE_TYPE] = Tags.HORIZONTAL_LAYER_STRUCTURE

all_structures["arbitrary_identifier"] = structure

simulation_settings[Tags.STRUCTURES] = all_structures

# ...
# Define further simulation settings
# ...

simulate(simulation_settings)

```

**create\_simulation\_volume ()** → dict

This method creates an in silico representation of a tissue as described in the settings file that is given.

**Parameters:** **settings** – a dictionary containing all relevant Tags for the simulation to be able to instantiate a tissue.

**Returns:** a path to a npz file containing characteristics of the simulated volume: absorption, scattering, anisotropy, oxygenation, and a segmentation mask. All of these are given as 3d numpy arrays.

### Segmentation-based volume creation

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```

class simpa.core.volume_creation_module.volume_creation_module_segmentation_based_adapter.
VolumeCreationModuleSegmentationBasedAdapter (global_settings:
simpa.utils.settings.Settings)

```

This volume creator expects a np.ndarray to be in the settings under the Tags.INPUT\_SEGMENTATION\_VOLUME tag and uses this array together with a SegmentationClass mapping which is a dict defined in the settings under Tags.SEGMENTATION\_CLASS\_MAPPING.

With this, an even greater utility is warranted.

**create\_simulation\_volume ()** → dict

This method will be called to create a simulation volume.

### Optical forward modelling

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```

class simpa.core.optical_simulation_module.OpticalForwardModuleBase (global_settings:
simpa.utils.settings.Settings)

```

Use this class as a base for implementations of optical forward models.

**abstract forward\_model** (absorption\_cm, scattering\_cm, anisotropy, illumination\_geometry, probe\_position\_mm)

A deriving class needs to implement this method according to its model.

**Parameters:**

- **absorption\_cm** – Absorption in units of per centimeter
- **scattering\_cm** – Scattering in units of per centimeter
- **anisotropy** – Dimensionless scattering anisotropy
- **illumination\_geometry** – A device that represents a detection geometry

**Returns:** Fluence in units of J/cm<sup>2</sup>

**run (device)**

Executes the respective simulation module

**Parameters:** **digital\_device\_twin** – The digital twin that can be used by the digital device\_twin.

## *mcx integration*

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.optical_simulation_module.optical_forward_model_mcx_adapter.OpticalForwardModelMcxAdapter(global_settings: simpa.utils.settings.Settings)
```

This class implements a bridge to the mcx framework to integrate mcx into SIMPA. MCX is a GPU-enabled Monte-Carlo model simulation of photon transport in tissue:

```
Fang, Qianqian, and David A. Boas. "Monte Carlo simulation of photon migration in 3D  
turbid media accelerated by graphics processing units."  
Optics express 17.22 (2009): 20178-20190.
```

```
forward_model (absorption_cm, scattering_cm, anisotropy, illumination_geometry,  
probe_position_mm)
```

A deriving class needs to implement this method according to its model.

**Parameters:**

- **absorption\_cm** – Absorption in units of per centimeter
- **scattering\_cm** – Scattering in units of per centimeter
- **anisotropy** – Dimensionless scattering anisotropy
- **illumination\_geometry** – A device that represents a detection geometry

**Returns:** Fluence in units of J/cm<sup>2</sup>

## *Acoustic forward modelling*

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.acoustic_forward_module.AcousticForwardModelBaseAdapter  
(global_settings: simpa.utils.settings.Settings)
```

This method is the entry method for running an acoustic forward model. It is invoked in the *simpa.core.simulation.simulate* method, but can also be called individually for the purposes of performing acoustic forward modeling only or in a different context.

The concrete will be chosen based on the:

```
Tags.ACOUSTIC_MODEL
```

tag in the settings dictionary.

**Parameters:** **settings** – The settings dictionary containing key-value pairs that determine the simulation. Here, it must contain the Tags.ACOUSTIC\_MODEL tag and any tags that might be required by the specific acoustic model.

**Raises:** **AssertionError** – an assertion error is raised if the Tags.ACOUSTIC\_MODEL tag is not given or points to an unknown acoustic forward model.

**abstract forward\_model** (detection\_geometry) → numpy.ndarray

This method performs the acoustic forward modeling given the initial pressure distribution and the acoustic tissue properties contained in the settings file. A deriving class needs to implement this method according to its model.

**Returns:** time series pressure data

**run** (digital\_device\_twin)

Call this method to invoke the simulation process.

**Parameters:** **digital\_device\_twin** –

**Returns:** a numpy array containing the time series pressure data per detection element

## *k-Wave integration*

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI)  
 SPDX-License-Identifier: MIT

```
class simpa.core.acoustic_forward_module.acoustic_forward_module_k_wave_adapter.AcousticForwardModelKWaveAdapter (global_settings: simpa.utils.settings.Settings)
```

The KwaveAcousticForwardModel adapter enables acoustic simulations to be run with the k-wave MATLAB toolbox. k-Wave is a free toolbox (<http://www.k-wave.org/>) developed by Bradley Treeby and Ben Cox (University College London) and Jiri Jaros (Brno University of Technology).

In order to use this toolbox, MATLAB needs to be installed on your system and the path to the MATLAB binary needs to be specified in the settings dictionary.

In order to use the toolbox from with SIMPA, a number of parameters have to be specified in the settings dictionary:

```
The initial pressure distribution:
    Tags.OPTICAL_MODEL_INITIAL_PRESSURE
Acoustic tissue properties:
    Tags.PROPERTY_SPEED_OF_SOUND
    Tags.PROPERTY_DENSITY
    Tags.PROPERTY_ALPHA_COEFF
The digital twin of the imaging device:
    Tags.DIGITAL_DEVICE
Other parameters:
    Tags.PERFORM_UPSAMPLING
    Tags.SPACING_MM
    Tags.UPSCALE_FACTOR
    Tags.PROPERTY_ALPHA_POWER
    Tags.GPU
    Tags.PMLInside
    Tags.PMLAlpha
    Tags.PlotPML
    Tags.RECORDMOVIE
    Tags.MOVIE_NAME
    Tags.ACOUSTIC_LOG_SCALE
    Tags.SENSOR_DIRECTIVITY_PATTERN
```

Many of these will be set automatically by SIMPA, but you may use the `simpa.utils.settings_generator` convenience methods to generate settings files that contain sensible defaults for these parameters.

Please also refer to the `simpa_examples` scripts to see how the settings file can be parametrized successfully.

**forward\_model** (detection\_geometry) → numpy.ndarray

This method performs the acoustic forward modeling given the initial pressure distribution and the acoustic tissue properties contained in the settings file. A deriving class needs to implement this method according to its model.

**Returns:** time series pressure data

## Image reconstruction

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

`class simpa.core.reconstruction_module.ReconstructionAdapterBase (global_settings: simpa.utils.settings.Settings)`

This class is the main entry point to perform image reconstruction using the SIMPA toolkit. All information necessary for the respective reconstruction method must be contained in the respective settings dictionary.

**abstract reconstruction\_algorithm** (time\_series\_sensor\_data, detection\_geometry: `simpa.core.device_digital_twins.devices.detection_geometries.detection_geometry_base.DetectionGeometryBase`) → `numpy.ndarray`

A deriving class needs to implement this method according to its model.

**Parameters:**

- **time\_series\_sensor\_data** – the time series sensor data
- **detection\_geometry** –

**Returns:** a reconstructed photoacoustic image

**run** (device)

Executes the respective simulation module

**Parameters:** **digital\_device\_twin** – The digital twin that can be used by the digital device\_twin.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

`simpa.core.reconstruction_module.reconstruction_utils.apply_b_mode` (data: `numpy.ndarray` = None, method: `str` = None) → `numpy.ndarray`

Applies B-Mode specified method to data. Method is either envelope detection using hilbert transform (Tags.RECONSTRUCTION\_BMODE\_METHOD\_HILBERT\_TRANSFORM), absolute value (Tags.RECONSTRUCTION\_BMODE\_METHOD\_ABS) or none if nothing is specified is performed.

**Parameters:**

- **data** – (numpy array) data used for applying B-Mode method
- **method** – (str)  
 Tags.RECONSTRUCTION\_BMODE\_METHOD\_HILBERT\_TRANSFORM or  
 Tags.RECONSTRUCTION\_BMODE\_METHOD\_ABS

**Returns:** (numpy array) data with B-Mode method applied, all

`simpa.core.reconstruction_module.reconstruction_utils.bandpass_filtering` (data: None, `_VariableFunctionsClass.tensor` = None, time\_spacing\_in\_ms: `float` = None, cutoff\_lowpass: `int` = 8000000, cutoff\_highpass: `int` = 100000, tukey\_alpha: `float` = 0.5) → `None._VariableFunctionsClass.tensor`

Apply a bandpass filter with cutoff values at *cutoff\_lowpass* and *cutoff\_highpass* MHz and a tukey window with alpha value of *tukey\_alpha* inbetween on the *data* in Fourier space.

**Parameters:**

- **data** – (torch tensor) data to be filtered
- **time\_spacing\_in\_ms** – (float) time spacing in milliseconds, e.g. 2.5e-5
- **cutoff\_lowpass** – (int) Signal above this value will be ignored (in MHz)
- **cutoff\_highpass** – (int) Signal below this value will be ignored (in MHz)
- **tukey\_alpha** – (float) transition value between 0 (rectangular) and 1 (Hann window)

**Returns:** (torch tensor) filtered data

`simpa.core.reconstruction_module.reconstruction_utils.get_apodization_factor` (apodization\_method: `str` = 'BoxApodization', dimensions: `tuple` = None, n\_sensor\_elements=None, device: `torch.device` = 'cpu') → `None._VariableFunctionsClass.tensor`  
 Construct apodization factors according to *apodization\_method* [hann, hamming or box apodization (default)] for given dimensions and *n\_sensor\_elements*.

**Parameters:**

- **apodization\_method** – (str) Apodization method, one of Tags.RECONSTRUCTION\_APODIZATION\_HANN, Tags.RECONSTRUCTION\_APODIZATION\_HAMMING and Tags.RECONSTRUCTION\_APODIZATION\_BOX (default)
- **dimensions** – (tuple) size of each dimension of reconstructed image as int, might have 2 or 3 entries.
- **n\_sensor\_elements** – (int) number of sensor elements
- **device** – (torch device) PyTorch tensor device

**Returns:** (torch tensor) tensor with apodization factors which can be multiplied with DAS values

```
simpa.core.reconstruction_module.reconstruction_utils.reconstruction_mode_transformation
(time_series_sensor_data: None._VariableFunctionsClass.tensor = None, mode: str =
'pressure') → None._VariableFunctionsClass.tensor
Transforms time_series_sensor_data for other modes, for example
Tags.RECONSTRUCTION_MODE_DIFFERENTIAL. Default mode is
Tags.RECONSTRUCTION_MODE_PRESSURE.
```

**Parameters:**

- **time\_series\_sensor\_data** – (torch tensor) Time series data to be transformed
- **mode** – (str) reconstruction mode: Tags.RECONSTRUCTION\_MODE\_PRESSURE (default) or Tags.RECONSTRUCTION\_MODE\_DIFFERENTIAL

**Returns:** (torch tensor) potentially transformed tensor

## Backprojection

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.reconstruction_module.reconstruction_module_delay_and_sum_adapter.Image
ReconstructionModuleDelayAndSumAdapter (global_settings:
simpa.utils.settings.Settings)
```

```
reconstruction_algorithm(time_series_sensor_data, detection_geometry: simpa.core.device_digital_twins.devices.detection_geometries.detection_geometry_base.DetectionGeometryBase)
```

Applies the Delay and Sum beamforming algorithm [1] to the time series sensor data (2D numpy array where the first dimension corresponds to the sensor elements and the second to the recorded time steps) with the given beamforming settings (dictionary). A reconstructed image (2D numpy array) is returned. This implementation uses PyTorch Tensors to perform computations and is able to run on GPUs.

[1] T. Kirchner et al. 2018, “Signed Real-Time Delay Multiply and Sum Beamforming for Multispectral Photoacoustic Imaging”, <https://doi.org/10.3390/jimaging4100121>

```
simpa.core.reconstruction_module.reconstruction_module_delay_and_sum_adapter.reconstruct_delay_and_sum
(time_series_sensor_data: numpy.ndarray, detection_geometry: simpa.core.device_digital_twins.devices.detection_geometries.detection_geometry_base.DetectionGeometryBase,
settings: dict = None, sound_of_speed: int = 1540, time_spacing: float = 2.5e-08,
sensor_spacing: float = 0.1) → numpy.ndarray
```

Convenience function for reconstructing time series data using Delay and Sum algorithm implemented in PyTorch

**Parameters:**

- **time\_series\_sensor\_data** – (2D numpy array) sensor data of shape (sensor elements, time steps)
- **detection\_geometry** – The DetectionGeometryBase that should be used to reconstruct the given time series data
- **settings** – (dict) settings dictionary: by default there is none and the other parameters are used instead, but if parameters are given in the settings those will be used instead of parsed arguments)
- **sound\_of\_speed** – (int) speed of sound in medium in meters per second (default: 1540 m/s)
- **time\_spacing** – (float) time between sampling points in seconds (default: 2.5e-8 s which is equal to 40 MHz)
- **sensor\_spacing** – (float) space between sensor elements in millimeters (default: 0.1 mm)

**Returns:** (2D numpy array) reconstructed image as 2D numpy array

### Delay-Multiply-And-Sum (DMAS)

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.reconstruction_module.reconstruction_module_delay_multiply_and_sum_adapter.ImageReconstructionModuleDelayMultiplyAndSumAdapter (global_settings:
simpa.utils.settings.Settings)
```

```
reconstruction_algorithm(time_series_sensor_data, detection_geometry: simpa.core.device_digital_twins.devices.detection_geometries.detection_geometry_base.DetectionGeometryBase)
```

Applies the Delay Multiply and Sum beamforming algorithm [1] to the time series sensor data (2D numpy array where the first dimension corresponds to the sensor elements and the second to the recorded time steps) with the given beamforming settings (dictionary). A reconstructed image (2D numpy array) is returned. This implementation uses PyTorch Tensors to perform computations and is able to run on GPUs.

[1] T. Kirchner et al. 2018, "Signed Real-Time Delay Multiply and Sum Beamforming for Multispectral Photoacoustic Imaging", <https://doi.org/10.3390/jimaging4100121>

```
simpa.core.reconstruction_module.reconstruction_module_delay_multiply_and_sum_adapter.reconstruction_algorithm(
time_series_sensor_data: numpy.ndarray, detection_geometry: simpa.core.device_digital_twins.devices.detection_geometries.detection_geometry_base.DetectionGeometryBase,
settings: dict = None, sound_of_speed: int = 1540, time_spacing: float = 2.5e-08,
sensor_spacing: float = 0.1) → numpy.ndarray
```

Convenience function for reconstructing time series data using Delay and Sum algorithm implemented in PyTorch

**Parameters:**

- **time\_series\_sensor\_data** – (2D numpy array) sensor data of shape (sensor elements, time steps)
- **detection\_geometry** – The DetectionGeometryBase to use for the reconstruction of the given time series data
- **settings** – (dict) settings dictionary: by default there is none and the other parameters are used instead, but if parameters are given in the settings those will be used instead of parsed arguments)
- **sound\_of\_speed** – (int) speed of sound in medium in meters per second (default: 1540 m/s)
- **time\_spacing** – (float) time between sampling points in seconds (default: 2.5e-8 s which is equal to 40 MHz)
- **sensor\_spacing** – (float) space between sensor elements in millimeters (default: 0.1 mm)



**Returns:** (2D numpy array) reconstructed image as 2D numpy array

### ***signed Delay-Multiply-And-Sum (sDMAS)***

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.reconstruction_module.reconstruction_module_signed_delay_multiply_and_sum_adapter.ImageReconstructionModuleSignedDelayMultiplyAndSumAdapter
(global_settings: simpa.utils.settings.Settings)
```

**reconstruction\_algorithm**(time\_series\_sensor\_data, detection\_geometry: simpa.core.device\_digital\_twins.devices.detection\_geometries.detection\_geometry\_base.DetectionGeometryBase)

Applies the signed Delay Multiply and Sum beamforming algorithm [1] to the time series sensor data (2D numpy array where the first dimension corresponds to the sensor elements and the second to the recorded time steps) with the given beamforming settings (dictionary). A reconstructed image (2D numpy array) is returned. This implementation uses PyTorch Tensors to perform computations and is able to run on GPUs.

[1] T. Kirchner et al. 2018, "Signed Real-Time Delay Multiply and Sum Beamforming for Multispectral Photoacoustic Imaging", <https://doi.org/10.3390/jimaging4100121>

```
simpa.core.reconstruction_module.reconstruction_module_signed_delay_multiply_and_sum_adapter.
(time_series_sensor_data: numpy.ndarray, detection_geometry: simpa.core.device_digital_twins.devices.detection_geometries.detection_geometry_base.DetectionGeometryBase,
settings: dict = None, sound_of_speed: int = 1540, time_spacing: float = 2.5e-08,
sensor_spacing: float = 0.1) → numpy.ndarray
```

Convenience function for reconstructing time series data using Delay and Sum algorithm implemented in PyTorch

#### **Parameters:**

- **time\_series\_sensor\_data** – (2D numpy array) sensor data of shape (sensor elements, time steps)
- **settings** – (dict) settings dictionary: by default there is none and the other parameters are used instead, but if parameters are given in the settings those will be used instead of parsed arguments)
- **sound\_of\_speed** – (int) speed of sound in medium in meters per second (default: 1540 m/s)
- **time\_spacing** – (float) time between sampling points in seconds (default: 2.5e-8 s which is equal to 40 MHz)
- **sensor\_spacing** – (float) space between sensor elements in millimeters (default: 0.1 mm)

**Returns:** (2D numpy array) reconstructed image as 2D numpy array

### ***Time Reversal***

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.reconstruction_module.reconstruction_module_time_reversal_adapter.ReconstructionModuleTimeReversalAdapter (global_settings: simpa.utils.settings.Settings)
```

The time reversal adapter includes the time reversal reconstruction algorithm implemented by the k-Wave toolkit into SIMPA.

Time reversal reconstruction uses the time series data and computes the forward simulation model backwards in time:

```
Treeby, Bradley E., Edward Z. Zhang, and Benjamin T. Cox.
"Photoacoustic tomography in absorbing acoustic media using
time reversal." Inverse Problems 26.11 (2010): 115003.
```

```
get_acoustic_properties(input_data: dict, detection_geometry)
```



This method extracts the acoustic tissue properties from the settings dictionary and amends the information to the input\_data.

**Parameters:**

- **global\_settings** – the settings dictionary containing key value pairs with the simulation instructions.
- **input\_data** – a dictionary containing the information needed for time reversal.

**reconstruction\_algorithm**(time\_series\_sensor\_data, detection\_geometry)

A deriving class needs to implement this method according to its model.

**Parameters:**

- **time\_series\_sensor\_data** – the time series sensor data
- **detection\_geometry** –

**Returns:** a reconstructed photoacoustic image

**reorder\_time\_series\_data**(time\_series\_sensor\_data, detection\_geometry)

Reorders the time series data to match the order that is assumed by kwave during image reconstruction with TimeReversal.

The main issue here is, that, while forward modelling allows for the definition of 3D cuboid bounding boxes for the detector elements, TimeReversal does not implement this feature. Instead, a binary mask is given and these are indexed in a column-row-wise manner in the output. The default np.argsort() method does not yield the same result as expected by k-Wave. Hence, this workaround.

## Processing Components

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.processing_components.ProcessingComponent (global_settings,
component_settings_key: str)
```

Defines a simulation component, which can be used to pre- or post-process simulation data.

## Noise Models

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.processing_components.noise.gamma_noise.GammaNoiseProcessingComponent
(global_settings, component_settings_key: str)
```

Applies Gaussian noise to the defined data field. The noise will be applied to all wavelengths.

Component Settings:

```
Tags.NOISE_SHAPE (default: 2)
Tags.NOISE_SCALE (default: 2)
Tags.NOISE_MODE (default: Tags.NOISE_MODE_ADDITIVE)
Tags.DATA_FIELD (required)
```

**run**(device)

Executes the respective simulation module

**Parameters:** **digital\_device\_twin** – The digital twin that can be used by the digital device\_twin.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.processing_components.noise.gaussian_noise.GaussianNoiseProcessingComponent
(global_settings, component_settings_key: str)
```

Applies Gaussian noise to the defined data field. The noise will be applied to all wavelengths. Component Settings:

```
Tags.NOISE_MEAN (default: 0)
Tags.NOISE_STD (default: 1)
Tags.NOISE_MODE (default: Tags.NOISE_MODE_ADDITIVE)
Tags.NOISE_NON_NEGATIVITY_CONSTRAINT (default: False)
Tags.DATA_FIELD (required)
```

**run (device)**

Executes the respective simulation module

**Parameters:** **digital\_device\_twin** – The digital twin that can be used by the digital device\_twin.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

*class*

`simpa.core.processing_components.noise.poisson_noise.PoissonNoiseProcessingComponent`  
(global\_settings, component\_settings\_key: str)

Applies Gaussian noise to the defined data field. The noise will be applied to all wavelengths.

Component Settings:

```
Tags.NOISE_MEAN (default: 3)
Tags.NOISE_MODE (default: Tags.NOISE_MODE_ADDITIVE)
Tags.DATA_FIELD (required)
```

**run (device)**

Executes the respective simulation module

**Parameters:** **digital\_device\_twin** – The digital twin that can be used by the digital device\_twin.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

*class* `simpa.core.processing_components.noise.salt_and_pepper_noise.SaltAndPepperNoiseProcessingComponent`  
(global\_settings, component\_settings\_key: str)

Applies salt and pepper noise to the defined data field. The noise will be applied to all wavelengths.

The noise will be 50% salt and 50% pepper noise, but both can be set to the same value using the NOISE\_MIN and NOISE\_MAX fields.

Component Settings:

```
Tags.NOISE_MIN (default: min(data_field))
Tags.NOISE_MAX (default: max(data_field))
Tags.NOISE_FREQUENCY (default: 0.01)
Tags.DATA_FIELD (required)
```

**run (device)**

Executes the respective simulation module

**Parameters:** **digital\_device\_twin** – The digital twin that can be used by the digital device\_twin.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

*class*

`simpa.core.processing_components.noise.uniform_noise.UniformNoiseProcessingComponent`  
(global\_settings, component\_settings\_key: str)

Applies uniform noise to the defined data field. The noise will be applied to all wavelengths.

The noise will be uniformly distributed between [min, max].

Component Settings:

```
Tags.NOISE_MIN (default: 0)
Tags.NOISE_MAX (default: 1)
```

```
Tags.NOISE_MODE (default: Tags.NOISE_MODE_ADDITIVE)
Tags.DATA_FIELD (required)
```

**run (device)**

Executes the respective simulation module

**Parameters:** **digital\_device\_twin** – The digital twin that can be used by the digital device\_twin.

## Digital device twins

At every step along the forward simulation, knowledge of the photoacoustic device that is used for the measurements is needed. This is important to reflect characteristic artefacts and challenges for the respective device.

To this end, we have included digital twins of commonly used devices into the SIMPA core. Additionally, we have included detection geometries and illumination geometries that can be used to create custom photoacoustic devices for simulation.

## Detection Geometries

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.device_digital_twins.devices.detection_geometries.detection_geometry_base.DetectionGeometryBase (number_detector_elements, detector_element_width_mm, detector_element_length_mm, center_frequency_hz, bandwidth_percent, sampling_frequency_mhz, probe_width_mm, device_position_mm: numpy.ndarray = None)
```

This class represents an illumination geometry

**abstract** **get\_detector\_element\_orientations** (global\_settings: simpa.utils.settings.Settings) → numpy.ndarray

This method yields a normalised orientation vector for each detection element. The length of this vector is the same as the one obtained via the position methods:

```
get_detector_element_positions_base_mm
get_detector_element_positions_accounting_for_device_position_mm
```

**Returns:** a numpy array that contains normalised orientation vectors for each detection element

**get\_detector\_element\_positions\_accounting\_for\_device\_position\_mm ()** → numpy.ndarray

Similar to:

```
get_detector_element_positions_base_mm
```

This method returns the absolute positions of the detection elements relative to the device position in the imaged volume, where the device position is defined by the following tag:

```
Tags.DIGITAL_DEVICE_POSITION
```

**Returns:** A numpy array containing the coordinates of the detection elements

**get\_detector\_element\_positions\_accounting\_for\_field\_of\_view ()** → numpy.ndarray

Similar to:

```
get_detector_element_positions_base_mm
```

This method returns the absolute positions of the detection elements relative to the device position in the imaged volume, where the device position is defined by the following tag:

```
Tags.DIGITAL_DEVICE_POSITION
```

**Returns:** A numpy array containing the coordinates of the detection elements

**abstract** `get_detector_element_positions_base_mm ()` → `numpy.ndarray`

Defines the abstract positions of the detection elements in an arbitrary coordinate system. Typically, the center of the field of view is defined as the origin.

To obtain the positions in an interpretable coordinate system, please use the other method:

```
get_detector_element_positions_accounting_for_device_position_mm
```

**Returns:** A numpy array containing the position vectors of the detection elements.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.device_digital_twins.devices.detection_geometries.curved_array.CurvedArrayDetectionGeometry (pitch_mm=0.5, radius_mm=40, number_detector_elements=256,
detector_element_width_mm=0.24, detector_element_length_mm=13,
center_frequency_hz=3960000.0, bandwidth_percent=55, sampling_frequency_mhz=40,
angular_origin_offset=3.141592653589793, device_position_mm=None)
```

This class represents a digital twin of an ultrasound detection device with a curved detection geometry

**check\_settings\_prerequisites** (`global_settings: simpa.utils.settings.Settings`) → `bool`

It might be that certain device geometries need a certain dimensionality of the simulated PAI volume, or that it required the existence of certain Tags in the global `global_settings`. To this end, a PAI device should use this method to inform the user about a mismatch of the desired device and throw a `ValueError` if that is the case.

**Raises:** **ValueError** – raises a value error if the prerequisites are not matched.

**Returns:** True if the prerequisites are met, False if they are not met, but no exception has been raised.

**get\_detector\_element\_orientations** (`global_settings: simpa.utils.settings.Settings`) → `numpy.ndarray`

This method yields a normalised orientation vector for each detection element. The length of this vector is the same as the one obtained via the position methods:

```
get_detector_element_positions_base_mm
get_detector_element_positions_accounting_for_device_position_mm
```

**Returns:** a numpy array that contains normalised orientation vectors for each detection element

**get\_detector\_element\_positions\_base\_mm ()** → `numpy.ndarray`

Defines the abstract positions of the detection elements in an arbitrary coordinate system. Typically, the center of the field of view is defined as the origin.

To obtain the positions in an interpretable coordinate system, please use the other method:

```
get_detector_element_positions_accounting_for_device_position_mm
```

**Returns:** A numpy array containing the position vectors of the detection elements.

**get\_field\_of\_view\_extent\_mm ()** → `numpy.ndarray`

Returns the field of view extent of this imaging device. It is defined as a numpy array of the shape [xs, xe, ys, ye, zs, ze], where x, y, and z denote the coordinate axes and s and e denote the start and end positions. These coordinates are defined relatively to the probe position.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.device_digital_twins.devices.detection_geometries.linear_array.LinearArrayDetectionGeometry (pitch_mm=0.5, number_detector_elements=100,
detector_element_width_mm=0.24, detector_element_length_mm=0.5,
center_frequency_hz=3960000.0, bandwidth_percent=55, sampling_frequency_mhz=40,
device_position_mm: numpy.ndarray = None)
```

This class represents a digital twin of an ultrasound detection device with a linear detection geometry.

**check\_settings\_prerequisites** (global\_settings: `simpa.utils.settings.Settings`) → bool

It might be that certain device geometries need a certain dimensionality of the simulated PAI volume, or that it required the existence of certain Tags in the global `global_settings`. To this end, a PAI device should use this method to inform the user about a mismatch of the desired device and throw a `ValueError` if that is the case.

**Raises:** **ValueError** – raises a value error if the prerequisites are not matched.

**Returns:** True if the prerequisites are met, False if they are not met, but no exception has been raised.

**get\_detector\_element\_orientations** (global\_settings: `simpa.utils.settings.Settings`) → `numpy.ndarray`

This method yields a normalised orientation vector for each detection element. The length of this vector is the same as the one obtained via the position methods:

```
get_detector_element_positions_base_mm
get_detector_element_positions_accounting_for_device_position_mm
```

**Returns:** a numpy array that contains normalised orientation vectors for each detection element

**get\_detector\_element\_positions\_base\_mm** () → `numpy.ndarray`

Defines the abstract positions of the detection elements in an arbitrary coordinate system. Typically, the center of the field of view is defined as the origin.

To obtain the positions in an interpretable coordinate system, please use the other method:

```
get_detector_element_positions_accounting_for_device_position_mm
```

**Returns:** A numpy array containing the position vectors of the detection elements.

**get\_field\_of\_view\_extent\_mm** () → `numpy.ndarray`

Returns the field of view extent of this imaging device. It is defined as a numpy array of the shape [xs, xe, ys, ye, zs, ze], where x, y, and z denote the coordinate axes and s and e denote the start and end positions. These coordinates are defined relatively to the probe position.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI)  
SPDX-License-Identifier: MIT

```
class simpa.core.device_digital_twins.devices.detection_geometries.planar_array.PlanarArrayDetectionGeometry
    (pitch_mm=0.5, number_detector_elements_x=100,
     number_detector_elements_y=100, detector_element_width_mm=0.24,
     detector_element_length_mm=0.5, center_frequency_hz=3960000.0, bandwidth_percent=55,
     sampling_frequency_mhz=40, device_position_mm: numpy.ndarray = None)
```

This class represents a digital twin of a ultrasound detection device with a linear detection geometry.

**check\_settings\_prerequisites** (global\_settings: `simpa.utils.settings.Settings`) → bool

It might be that certain device geometries need a certain dimensionality of the simulated PAI volume, or that it required the existence of certain Tags in the global `global_settings`. To this end, a PAI device should use this method to inform the user about a mismatch of the desired device and throw a `ValueError` if that is the case.

**Raises:** **ValueError** – raises a value error if the prerequisites are not matched.

**Returns:** True if the prerequisites are met, False if they are not met, but no exception has been raised.

**get\_detector\_element\_orientations** (global\_settings: `simpa.utils.settings.Settings`) → `numpy.ndarray`

This method yields a normalised orientation vector for each detection element. The length of this vector is the same as the one obtained via the position methods:

```
get_detector_element_positions_base_mm
get_detector_element_positions_accounting_for_device_position_mm
```

**Returns:** a numpy array that contains normalised orientation vectors for each detection element

**get\_detector\_element\_positions\_base\_mm ()** → numpy.ndarray

Defines the abstract positions of the detection elements in an arbitrary coordinate system. Typically, the center of the field of view is defined as the origin.

To obtain the positions in an interpretable coordinate system, please use the other method:

```
get_detector_element_positions_accounting_for_device_position_mm
```

**Returns:** A numpy array containing the position vectors of the detection elements.

**get\_field\_of\_view\_extent\_mm ()** → numpy.ndarray

Returns the field of view extent of this imaging device. It is defined as a numpy array of the shape [xs, xe, ys, ye, zs, ze], where x, y, and z denote the coordinate axes and s and e denote the start and end positions. These coordinates are defined relatively to the probe position.

## Illumination Geometries

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.device_digital_twins.devices.illumination_geometries.illumination_geometry_base.IlluminationGeometryBase
```

This class represents an illumination geometry.

**check\_settings\_prerequisites (global\_settings: simpa.utils.settings.Settings)** → bool

It might be that certain device geometries need a certain dimensionality of the simulated PAI volume, or that it required the existence of certain Tags in the global global\_settings. To this end, a PAI device should use this method to inform the user about a mismatch of the desired device and throw a ValueError if that is the case.

**Raises:** **ValueError** – raises a value error if the prerequisites are not matched.

**Returns:** True if the prerequisites are met, False if they are not met, but no exception has been raised.

**abstract** **get\_mcx\_illuminator\_definition** (global\_settings: simpa.utils.settings.Settings, probe\_position\_mm: numpy.ndarray) → dict

IMPORTANT: This method creates a dictionary that contains tags as they are expected for the mcx simulation tool to represent the illumination geometry of this device.

**Parameters:**

- **global\_settings** – The global\_settings instance containing the simulation instructions
- **probe\_position\_mm** – the position of the probe in the volume

**Returns:** Dictionary that includes all parameters needed for mcx.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.device_digital_twins.devices.illumination_geometries.pencil_array_illumination.PencilArrayIlluminationGeometry (pitch_mm=0.5, number_illuminators_x=100, number_illuminators_y=100)
```

This class represents a slit illumination geometry. The device position is defined as the middle of the slit.

**get\_field\_of\_view\_extent\_mm ()** → numpy.ndarray

Returns the field of view extent of this imaging device. It is defined as a numpy array of the shape [xs, xe, ys, ye, zs, ze], where x, y, and z denote the coordinate axes and s and e denote the start and end positions. These coordinates are defined relatively to the probe position.

**get\_mcx\_illuminator\_definition** (global\_settings: simpa.utils.settings.Settings, probe\_position\_mm) → dict

IMPORTANT: This method creates a dictionary that contains tags as they are expected for the mcx simulation tool to represent the illumination geometry of this device.



**Parameters:**

- **global\_settings** – The `global_settings` instance containing the simulation instructions
- **probe\_position\_mm** – the position of the probe in the volume

**Returns:** Dictionary that includes all parameters needed for mcx.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.device_digital_twins.devices.illumination_geometries.pencil_beam_illumination.PencilBeamIlluminationGeometry
```

This class represents a pencil beam illumination geometry. The device position is defined as the middle of the slit.

**get\_field\_of\_view\_extent\_mm ()** → `numpy.ndarray`

Returns the field of view extent of this imaging device. It is defined as a numpy array of the shape [xs, xe, ys, ye, zs, ze], where x, y, and z denote the coordinate axes and s and e denote the start and end positions. These coordinates are defined relatively to the probe position.

**get\_mcx\_illuminator\_definition** (`global_settings:` `simpa.utils.settings.Settings`, `probe_position_mm`) → `dict`

IMPORTANT: This method creates a dictionary that contains tags as they are expected for the mcx simulation tool to represent the illumination geometry of this device.

**Parameters:**

- **global\_settings** – The `global_settings` instance containing the simulation instructions
- **probe\_position\_mm** – the position of the probe in the volume

**Returns:** Dictionary that includes all parameters needed for mcx.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.core.device_digital_twins.devices.illumination_geometries.slit_illumination.SlitIlluminationGeometry (slit_vector_mm: list = None, direction_vector_mm: list = None)
```

This class represents a slit illumination geometry. The device position is defined as the middle of the slit.

**get\_field\_of\_view\_extent\_mm ()** → `numpy.ndarray`

Returns the field of view extent of this imaging device. It is defined as a numpy array of the shape [xs, xe, ys, ye, zs, ze], where x, y, and z denote the coordinate axes and s and e denote the start and end positions. These coordinates are defined relatively to the probe position.

**get\_mcx\_illuminator\_definition** (`global_settings:` `simpa.utils.settings.Settings`, `probe_position_mm`) → `dict`

IMPORTANT: This method creates a dictionary that contains tags as they are expected for the mcx simulation tool to represent the illumination geometry of this device.

**Parameters:**

- **global\_settings** – The `global_settings` instance containing the simulation instructions
- **probe\_position\_mm** – the position of the probe in the volume

**Returns:** Dictionary that includes all parameters needed for mcx.

## Models of real world devices

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class
```

```
simpa.core.device_digital_twins.devices.pa_devices.ithera_msot_acuity.MSOTAcuityEcho (device_position_mm: numpy.ndarray = None)
```

This class represents a digital twin of the MSOT Acuity Echo, manufactured by iThera Medical, Munich, Germany (<https://www.ithera-medical.com/products/msot-acuity/>). It is based on the real specifications of the device, but due to the limitations of the possibilities how to represent a device in the software frameworks, constitutes only an approximation.

Some important publications that showcase the use cases of the MSOT Acuity and Acuity Echo device are:

Regensburger, Adrian P., et al. "Detection of collagens by multispectral optoacoustic tomography as an imaging biomarker for Duchenne muscular dystrophy." *Nature Medicine* 25.12 (2019): 1905-1915.

Knieling, Ferdinand, et al. "Multispectral Optoacoustic Tomography for Assessment of Crohn's Disease Activity." *The New England journal of medicine* 376.13 (2017): 1292.

**update\_settings\_for\_use\_of\_model\_based\_volume\_creator** (global\_settings: `simpa.utils.settings.Settings`)

This method can be overwritten by a photoacoustic device if the device poses special constraints to the volume that should be considered by the model-based volume creator.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

**class**

`simpa.core.device_digital_twins.devices.pa_devices.ithera_msot_invision.InVision256TF`  
(`device_position_mm: numpy.ndarray = None`)

This class represents a digital twin of the InVision 256-TF, manufactured by iThera Medical, Munich, Germany (<https://www.ithera-medical.com/products/msot-invision/>). It is based on the real specifications of the device, but due to the limitations of the possibilities how to represent a device in the software frameworks, constitutes only an approximation.

Some important publications that showcase the use cases of the InVision series devices are:

Joseph, James, et al. "Evaluation of precision in optoacoustic tomography for preclinical imaging in living subjects." *Journal of Nuclear Medicine* 58.5 (2017): 807-814.

Merkle, Elena, et al. "Whole-body live mouse imaging by hybrid reflection-mode ultrasound and optoacoustic tomography." *Optics letters* 40.20 (2015): 4643-4646.

**get\_field\_of\_view\_extent\_mm()** → `numpy.ndarray`

Returns the field of view extent of this imaging device. It is defined as a numpy array of the shape [xs, xe, ys, ye, zs, ze], where x, y, and z denote the coordinate axes and s and e denote the start and end positions. These coordinates are defined relatively to the probe position.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

**class** `simpa.core.device_digital_twins.devices.pa_devices.ithera_rsom.RSOMExplorerP50`  
(`element_spacing_mm=0.02, number_elements_x=10, number_elements_y=10, device_position_mm: numpy.ndarray = None`)

This class represents an approximation of the Raster-scanning Optoacoustic Mesoscopy (RSOM) device built by iThera Medical (Munich, Germany). Please refer to the company's website for more information (<https://www.ithera-medical.com/products/rsom-explorer-p50/>).

Since simulating thousands of individual forward modeling steps to obtain a single raster-scanned image is computationally not feasible, we approximate the process with a device design that has detection elements across the entire field of view. Because of this limitation we also need to approximate the light source with a homogeneous illumination across the field of view.

The digital device is modeled based on the reported specifications of the RSOM Explorer P50 system. Technical details of the system can be found in the dissertation of Mathias Schwarz (<https://mediatum.ub.tum.de/doc/1324031/1324031.pdf>) and you can find more details on use cases of the device in the following literature sources:



Yew, Yik Weng, et al. "Raster-scanning optoacoustic mesoscopy (RSOM) imaging as an objective disease severity tool in atopic dermatitis patients." *Journal of the American Academy of Dermatology* (2020).

Hindelang, B., et al. "Non-invasive imaging in dermatology and the unique potential of raster-scan optoacoustic mesoscopy." *Journal of the European Academy of Dermatology and Venereology* 33.6 (2019): 1051-1061.

## Module: utils

The utils module contains several general-purpose utility functions whose purpose it is to facilitate the use of SIMPA. The most important of these is the Tags class, which defines the strings and data types that have to be used for the keys and values of the settings dictionary.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

`simpa.utils.calculate.calculate_gruneisen_parameter_from_temperature`  
 (temperature\_in\_celcius)

This function returns the dimensionless gruneisen parameter based on a heuristic formula that was determined experimentally:

```
@book{wang2012biomedical,
  title={Biomedical optics: principles and imaging},
  author={Wang, Lihong V and Wu, Hsin-i},
  year={2012},
  publisher={John Wiley & Sons}
}
```

**Parameters:** `temperature_in_celcius` – the temperature in degrees celcius

**Returns:** a floating point number, if `temperature_in_celcius` is a number or a float array, if `temperature_in_celcius` is an array

`simpa.utils.calculate.calculate_oxygenation` (molecule\_list)

**Returns:** an oxygenation value between 0 and 1 if possible, or None, if not computable.

`simpa.utils.calculate.create_spline_for_range` (xmin\_mm=0, xmax\_mm=10,  
 maximum\_y\_elevation\_mm=1, spacing=0.1)

Creates a functional that simulates distortion along the y position between the minimum and maximum x positions. The elevation can never be smaller than 0 or bigger than `maximum_y_elevation_mm`.

**Parameters:**

- `xmin_mm` – the minimum x axis value the return functional is defined in
- `xmax_mm` – the maximum x axis value the return functional is defined in
- `maximum_y_elevation_mm` – the maximum y axis value the return functional will yield

**Returns:** a functional that describes a distortion field along the y axis

`simpa.utils.calculate.min_max_normalization` (data: `numpy.ndarray` = None) → `numpy.ndarray`  
 Normalizes the given data by applying min max normalization. The resulting array has values between 0 and 1 inclusive.

**Parameters:** `data` – (numpy array) data to be normalized

**Returns:** (numpy array) normalized array

`simpa.utils.calculate.randomize_uniform` (min\_value: float, max\_value: float)  
 returns a uniformly drawn random number in [min\_value, max\_value[

**Parameters:**

- **min\_value** – minimum value
- **max\_value** – maximum value

**Returns:** random number in [min\_value, max\_value[

`simpa.utils.calculate.rotation` (angles)

Rotation matrix around the x-, y-, and z-axis with angles [theta\_x, theta\_y, theta\_z].

**Parameters:** **angles** – Angles through which the matrix is supposed to rotate in the form of [theta\_x, theta\_y, theta\_z].

**Returns:** rotation matrix

`simpa.utils.calculate.rotation_matrix_between_vectors` (a, b)

Returns the rotation matrix from a to b

**Parameters:**

- **a** – 3D vector to rotate
- **b** – 3D target vector

**Returns:** rotation matrix

`simpa.utils.calculate.rotation_x` (theta)

Rotation matrix around the x-axis with angle theta.

**Parameters:** **theta** – Angle through which the matrix is supposed to rotate.

**Returns:** rotation matrix

`simpa.utils.calculate.rotation_y` (theta)

Rotation matrix around the y-axis with angle theta.

**Parameters:** **theta** – Angle through which the matrix is supposed to rotate.

**Returns:** rotation matrix

`simpa.utils.calculate.rotation_z` (theta)

Rotation matrix around the z-axis with angle theta.

**Parameters:** **theta** – Angle through which the matrix is supposed to rotate.

**Returns:** rotation matrix

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

`simpa.utils.constants.EPS = 1e-20`

Defines the smallest increment that should be considered by SIMPA.

`class simpa.utils.constants.SaveFilePaths`

The save file paths specify the path of a specific data structure in the dictionary of the simpa output hdf5. All of these paths have to be used like: SaveFilePaths.PATH + "data\_structure"

`class simpa.utils.constants.SegmentationClasses`

The segmentation classes define which "tissue types" are modelled in the simulation volumes.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

`simpa.utils.deformation_manager.create_deformation_settings` (bounds\_mm,  
maximum\_z\_elevation\_mm=1, filter\_sigma=1, cosine\_scaling\_factor=4)

FIXME

`simpa.utils.deformation_manager.get_functional_from_deformation_settings`

(deformation\_settings: dict)

FIXME

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
simpa.utils.dict_path_manager.generate_dict_path(data_field, wavelength: (<class 'int'>, <class 'float'>) = None) → str
```

Generates a path within an hdf5 file in the SIMPA convention

**Parameters:**

- **data\_field** – Data field that is supposed to be stored in an hdf5 file.
- **wavelength** – Wavelength of the current simulation.

**Returns:** String which defines the path to the data\_field.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI)  
SPDX-License-Identifier: MIT

```
class simpa.utils.path_manager.PathManager(environment_path=None)
```

As a pipelining tool that serves as a communication layer between different numerical forward models and processing tools, SIMPA needs to be configured with the paths to these tools on your local hard drive. To this end, we have implemented the *PathManager* class that you can import to your project using *from simpa.utils import PathManager*. The *PathManager* looks for a *path\_config.env* file (just like the one we provided in the *simpa\_examples*) in the following places in this order:

1. The optional path you give the *PathManager*
2. Your \$HOME\$ directory
3. The current working directory
4. The SIMPA home directory path

```
detect_local_path_config ()
```

This methods looks in the default local paths for a *path\_config.env* file.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI)  
SPDX-License-Identifier: MIT

```
class simpa.utils.settings.Settings(dictionary: dict = None)
```

The *Settings* class is a dictionary that contains all relevant settings for running a simulation in the SIMPA toolkit. It includes an automatic sanity check for input parameters using the *simpa.utils.Tags* class.

Usage: *Seetings*({Tags.KEY1: value1, Tags.KEY2: value2, ...})

```
get_acoustic_settings ()
```

" Returns the settings for the acoustic forward model that are saved in this settings dictionary

```
get_optical_settings ()
```

" Returns the settings for the optical forward model that are saved in this settings dictionary

```
get_reconstruction_settings ()
```

" Returns the settings for the reconstruction model that are saved in this settings dictionary

```
get_volume_creation_settings ()
```

" Returns the settings for the optical forward model that are saved in this settings dictionary

```
set_acoustic_settings (acoustic_settings: dict)
```

Replaces the currently stored acoustic forward model settings with the given dictionary

**Parameters:** **acoustic\_settings** – a dictionary containing the acoustic model settings

```
set_optical_settings (optical_settings: dict)
```

Replaces the currently stored optical settings with the given dictionary

**Parameters:** **optical\_settings** – a dictionary containing the optical settings

```
set_reconstruction_settings (reconstruction_settings: dict)
```

Replaces the currently stored reconstruction model settings with the given dictionary

**Parameters:** **reconstruction\_settings** – a dictionary containing the reconstruction model settings

**set\_volume\_creation\_settings**(volume\_settings: dict)

Replaces the currently stored volume creation settings with the given dictionary

**Parameters:** **volume\_settings** – a dictionary containing the volume creator settings

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI)  
SPDX-License-Identifier: MIT

**class simpa.utils.tags.Tags**

This class contains all 'Tags' for the use in the settings dictionary as well as strings that are used in SIMPA as naming conventions. Every Tag that is intended to be used as a key in the settings dictionary is represented by a tuple. The first element of the tuple is a string that corresponds to the name of the Tag. The second element of the tuple is a data type or a tuple of data types. The values that are assigned to the keys in the settings should match these data types. Their usage within the SIMPA package is divided in "SIMPA package", "module X", "adapter Y", "class Z" and "naming convention".

**ACOUSTIC\_LOG\_SCALE** = ('acoustic\_log\_scale', (<class 'bool'>, <class 'bool'>, <class 'numpy.bool\_'>))

If True, the movie of the kwave simulation will be recorded in a log scale.

Usage: adapter KwaveAcousticForwardModel

**ACOUSTIC\_MODEL** = ('acoustic\_model', <class 'str'>)

Choice of the used acoustic model.

Usage: module acoustic\_forward\_module

**ACOUSTIC\_MODEL\_BINARY\_PATH** = ('acoustic\_model\_binary\_path', <class 'str'>)

Absolute path of the location of the acoustic forward model binary.

Usage: module optical\_simulation\_module

**ACOUSTIC\_MODEL\_K\_WAVE** = 'kwave'

Corresponds to the kwave simulation.

Usage: module acoustic\_forward\_module, naming convention

**ACOUSTIC\_MODEL\_OUTPUT\_NAME** = 'acoustic\_forward\_model\_output'

Name of the acoustic forward model output field in the SIMPA output file.

Usage: naming convention

**ACOUSTIC\_MODEL\_SETTINGS** = ('acoustic\_model\_settings', <class 'dict'>)

Acoustic model settings.

**ACOUSTIC\_MODEL\_TEST** = 'simpa\_tests'

Corresponds to an adapter for testing purposes only.

Usage: module acoustic\_forward\_module, naming convention

**ACOUSTIC\_SIMULATION\_3D** = ('acoustic\_simulation\_3d', <class 'bool'>)

If True, simulates the acoustic forward model in 3D.

Usage: SIMPA package

**ADHERE\_TO\_DEFORMATION** = ('adhere\_to\_deformation', <class 'bool'>)

If True, a structure will be shifted according to the deformation.

Usage: adapter versatile\_volume\_creation

**BACKGROUND** = 'Background'

Corresponds to the name of a structure.

Usage: adapter versatile\_volume\_creation, naming convention

**BANDPASS\_CUTOFF\_HIGHPASS** = ('bandpass\_cutoff\_highpass', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Sets the cutoff threshold in MHz for highpass filtering, i.e. lower limit of the tukey filter. Default is 0.1 MHz

Usage: adapter PyTorchDASAdapter

**BANDPASS\_CUTOFF\_LOWPASS** = ('bandpass\_cutoff\_lowpass', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))  
Sets the cutoff threshold in MHz for lowpass filtering, i.e. upper limit of the tukey filter. Default is 8 MHz  
Usage: adapter PyTorchDASAdapter

**CIRCULAR\_TUBULAR\_STRUCTURE** = 'CircularTubularStructure'  
Corresponds to the CircularTubularStructure in the structure\_library.  
Usage: module volume\_creation\_module, naming\_convention

**CONSIDER\_PARTIAL\_VOLUME** = ('consider\_partial\_volume', <class 'bool'>)  
If True, the structure will be generated with its edges only occupying a partial volume of the voxel.  
Usage: adapter versatile\_volume\_creation

**DATA\_FIELD** = 'data\_field'  
Defines which data field a certain function shall be applied to.  
Usage: module core.processing\_components

**DEFORMATION\_X\_COORDINATES\_MM** = 'deformation\_x\_coordinates'  
Mesh that defines the x coordinates of the deformation.  
Usage: adapter versatile\_volume\_creation, naming\_convention

**DEFORMATION\_Y\_COORDINATES\_MM** = 'deformation\_y\_coordinates'  
Mesh that defines the y coordinates of the deformation.  
Usage: adapter versatile\_volume\_creation, naming\_convention

**DEFORMATION\_Z\_ELEVATIONS\_MM** = 'deformation\_z\_elevation'  
Mesh that defines the z coordinates of the deformation.  
Usage: adapter versatile\_volume\_creation, naming\_convention

**DEFORMED\_LAYERS\_SETTINGS** = ('deformed\_layers\_settings', <class 'dict'>)  
Settings that contain the functional which defines the deformation of the layers.  
Usage: adapter versatile\_volume\_creation

**DETECTOR\_ELEMENT\_WIDTH\_MM** = 'detector\_element\_width\_mm'  
Width of a detector element. Corresponds to the pitch - the distance between two detector element borders.  
Usage: module acoustic\_forward\_module, naming\_convention

**DIGITAL\_DEVICE** = ('digital\_device', <class 'str'>)  
Digital device that is chosen as illumination source and detector for the simulation.  
Usage: SIMPA package

**DIGITAL\_DEVICE\_MSOT\_ACUITY** = 'digital\_device\_msot'  
Corresponds to the MSOTAcuityEcho device.  
Usage: SIMPA package, naming\_convention

**DIGITAL\_DEVICE\_MSOT\_INVISION** = 'digital\_device\_invision'  
Corresponds to the InVision 256-TF device.  
Usage: SIMPA package, naming\_convention

**DIGITAL\_DEVICE\_POSITION** = ('digital\_device\_position', (<class 'list'>, <class 'tuple'>, <class 'numpy.ndarray'>))  
Position in [x, y, z] coordinates of the device in the generated volume.  
Usage: SIMPA package

**DIGITAL\_DEVICE\_RSOM** = 'digital\_device\_rsom'  
Corresponds to the RSOMExplorerP50 device.  
Usage: SIMPA package, naming\_convention

**DIGITAL\_DEVICE\_SLIT\_ILLUMINATION\_LINEAR\_DETECTOR** = 'digital\_device\_slit\_illumination\_linear\_detector'  
Corresponds to a PA device with a slit as illumination and a linear array as detection geometry.

Usage: SIMPA package, naming convention

**DIM\_VOLUME\_X\_MM** = ('volume\_x\_dim\_mm', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Extent of the x-axis of the generated volume.

Usage: SIMPA package

**DIM\_VOLUME\_Y\_MM** = ('volume\_y\_dim\_mm', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Extent of the y-axis of the generated volume.

Usage: SIMPA package

**DIM\_VOLUME\_Z\_MM** = ('volume\_z\_dim\_mm', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Extent of the z-axis of the generated volume.

Usage: SIMPA package

**DOWNSCALE\_FACTOR** = ('downscale\_factor', (<class 'int'>, <class 'float'>, <class 'numpy.int32'>, <class 'numpy.float64'>))

Downscale factor of the resampling in the qPAI reconstruction

Usage: module algorithms (iterative\_qPAI\_algorithm.py)

**ELLIPTICAL\_TUBULAR\_STRUCTURE** = 'EllipticalTubularStructure'

Corresponds to the EllipticalTubularStructure in the structure\_library.

Usage: module volume\_creation\_module, naming\_convention

**GPU** = ('gpu', (<class 'bool'>, <class 'bool'>, <class 'numpy.bool\_'>))

If True, uses all available gpu options of the used modules.

Usage: SIMPA package

**HORIZONTAL\_LAYER\_STRUCTURE** = 'HorizontalLayerStructure'

Corresponds to the HorizontalLayerStructure in the structure\_library.

Usage: module volume\_creation\_module, naming\_convention

**ILLUMINATION\_DIRECTION** = ('illumination\_direction', (<class 'list'>, <class 'tuple'>, <class 'numpy.ndarray'>))

Direction of the photon source as [x, y, z] vector used in mcx.

Usage: module optical\_modelling, adapter mcx\_adapter

**ILLUMINATION\_PARAM1** = ('illumination\_param1', (<class 'list'>, <class 'tuple'>, <class 'numpy.ndarray'>))

First parameter group of the specified illumination type as [x, y, z, w] vector used in mcx.

Usage: module optical\_modelling, adapter mcx\_adapter

**ILLUMINATION\_PARAM2** = ('illumination\_param2', (<class 'list'>, <class 'tuple'>, <class 'numpy.ndarray'>))

Second parameter group of the specified illumination type as [x, y, z, w] vector used in mcx.

Usage: module optical\_modelling, adapter mcx\_adapter

**ILLUMINATION\_POSITION** = ('illumination\_position', (<class 'list'>, <class 'tuple'>, <class 'numpy.ndarray'>))

Position of the photon source in [x, y, z] coordinates used in mcx.

Usage: module optical\_modelling, adapter mcx\_adapter

**ILLUMINATION\_TYPE** = ('optical\_model\_illumination\_type', <class 'str'>)

Type of the illumination geometry used in mcx.

Usage: module optical\_modelling, adapter mcx\_adapter

**ILLUMINATION\_TYPE\_DISK** = 'disk'

Corresponds to disk source in mcx.

Usage: adapter mcx\_adapter, naming\_convention

**ILLUMINATION\_TYPE\_DKFZ\_PAUS** = 'pasetup'

Corresponds to pasetup source in mcx. The geometrical definition is described in:

Usage: adapter mcx\_adapter, naming\_convention

**ILLUMINATION\_TYPE\_FOURIER** = 'fourier'

Corresponds to fourier source in mcx.

Usage: adapter mcx\_adapter, naming convention

**ILLUMINATION\_TYPE\_FOURIER\_X** = *'fourierx'*

Corresponds to fourierx source in mcx.

Usage: adapter mcx\_adapter, naming convention

**ILLUMINATION\_TYPE\_FOURIER\_X\_2D** = *'fourierx2d'*

Corresponds to fourierx2d source in mcx.

Usage: adapter mcx\_adapter, naming convention

**ILLUMINATION\_TYPE\_GAUSSIAN** = *'gaussian'*

Corresponds to gaussian source in mcx.

Usage: adapter mcx\_adapter, naming convention

**ILLUMINATION\_TYPE\_IPASC\_DEFINITION** = *'ipasc'*

Corresponds to a source definition in mcx.

Usage: adapter mcx\_adapter, naming convention

**ILLUMINATION\_TYPE\_MSOT\_ACUITY\_ECHO** = *'msot\_acuity\_echo'*

Corresponds to msot\_acuity\_echo source in mcx. The device is manufactured by iThera Medical, Munich, Germany (<https://www.ithera-medical.com/products/msot-acuity/>).

Usage: adapter mcx\_adapter, naming convention

**ILLUMINATION\_TYPE\_MSOT\_INVISION** = *'invision'*

Corresponds to a source definition in mcx.

Usage: adapter mcx\_adapter, naming convention

**ILLUMINATION\_TYPE\_PATTERN** = *'pattern'*

Corresponds to pattern source in mcx.

Usage: adapter mcx\_adapter, naming convention

**ILLUMINATION\_TYPE\_PATTERN\_3D** = *'pattern3d'*

Corresponds to pattern3d source in mcx.

Usage: adapter mcx\_adapter, naming convention

**ILLUMINATION\_TYPE\_PENCIL** = *'pencil'*

Corresponds to pencil source in mcx.

Usage: adapter mcx\_adapter, naming convention

**ILLUMINATION\_TYPE\_PENCILARRAY** = *'pencilarray'*

Corresponds to pencilarray source in mcx.

Usage: adapter mcx\_adapter, naming convention

**ILLUMINATION\_TYPE\_PLANAR** = *'planar'*

Corresponds to planar source in mcx.

Usage: adapter mcx\_adapter, naming convention

**ILLUMINATION\_TYPE\_RING** = *'ring'*

Corresponds to ring source in mcx.

Usage: adapter mcx\_adapter, naming convention

**ILLUMINATION\_TYPE\_SLIT** = *'slit'*

Corresponds to slit source in mcx.

Usage: adapter mcx\_adapter, naming convention

**IMAGE\_PROCESSING** = *'image\_processing'*

Location of the image algorithms outputs in the SIMPA output file.

Usage: naming convention

**INPUT\_SEGMENTATION\_VOLUME** = (*'input\_segmentation\_volume'*, <class 'numpy.ndarray'>)

Array that defines a segmented volume.

Usage: adapter segmentation\_based\_volume\_creator

**ITERATIVE\_RECONSTRUCTION\_CONSTANT\_REGULARIZATION** = ('constant\_regularization', (<class 'bool'>, <class 'bool'>, <class 'numpy.bool\_'>))

If True, the fluence regularization will be constant.

Usage: module algorithms (iterative\_qPAI\_algorithm.py)

**ITERATIVE\_RECONSTRUCTION\_MAX\_ITERATION\_NUMBER** = ('maximum\_iteration\_number', (<class 'int'>, <class 'numpy.integer'>))

Maximum number of iterations performed in iterative reconstruction if stopping criterion is not reached.

Usage: module algorithms (iterative\_qPAI\_algorithm.py)

**ITERATIVE\_RECONSTRUCTION\_REGULARIZATION\_SIGMA** = ('regularization\_sigma', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Sigma value used for constant regularization of fluence.

Usage: module algorithms (iterative\_qPAI\_algorithm.py)

**ITERATIVE\_RECONSTRUCTION\_SAVE\_INTERMEDIATE\_RESULTS** = ('save\_intermediate\_results', (<class 'bool'>, <class 'bool'>, <class 'numpy.bool\_'>))

If True, a list of all intermediate absorption updates (middle slices only) will be saved in a numpy file.

Usage: module algorithms (iterative\_qPAI\_algorithm.py)

**ITERATIVE\_RECONSTRUCTION\_STOPPING\_LEVEL** = ('iteration\_stopping\_level', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Ratio of improvement and preceding error at which iteration method stops. Usage: module algorithms (iterative\_qPAI\_algorithm.py)

**ITERATIVE\_qPAI\_RESULT** = 'iterative\_qpai\_result'

Name of the data field in which the iterative qPAI result will be stored.

Usage: naming convention

**K\_WAVE\_SPECIFIC\_DT** = ('dt\_acoustic\_sim', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Temporal resolution of kwave.

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter

**K\_WAVE\_SPECIFIC\_NT** = ('Nt\_acoustic\_sim', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Total time steps simulated by kwave.

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter

**LASER\_PULSE\_ENERGY\_IN\_MILLIJOULE** = ('laser\_pulse\_energy\_in\_millijoule', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>, <class 'list'>, <class 'range'>, <class 'tuple'>, <class 'numpy.ndarray'>))

Laser pulse energy used in the optical simulation.

Usage: module optical\_simulation\_module

**LINEAR\_UNMIXING\_RESULT** = 'linear\_unmixing\_result'

Name of the data field in which the linear unmixing result will be stored.

Usage: naming convention

**LOAD\_AND\_SAVE\_HDF5\_FILE\_AT\_THE\_END\_OF\_SIMULATION\_TO\_MINIMISE\_FILESIZE** = ('minimize\_file\_size', (<class 'bool'>, <class 'bool'>, <class 'numpy.bool\_'>))

If not set to False, the HDF5 file will be optimised after the simulations are done. Usage: simpa.core.simulation.simulate

**MAX\_DEFORMATION\_MM** = 'max\_deformation'

Maximum deformation in z-direction.

Usage: adapter versatile\_volume\_creation, naming convention

**MCX\_ASSUMED\_ANISOTROPY** = ('mcx\_seed', (<class 'float'>, <class 'float'>))



The anisotropy that should be assumed for the mcx simulations. If not set, a default value of 0.9 will be assumed. Usage: module optical\_modelling, adapter mcx\_adapter

**MCX\_SEED** = ('mcx\_seed', (<class 'int'>, <class 'numpy.integer'>))

Specific seed for random initialisation in mcx.

if not set, Tags.RANDOM\_SEED will be used instead. Usage: module optical\_modelling, adapter mcx\_adapter

**MEDIUM\_TEMPERATURE\_CELCIUS** = ('medium\_temperature', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Temperature of the simulated volume.

Usage: module noise\_simulation

**MOLECULE\_COMPOSITION** = ('molecule\_composition', <class 'list'>)

List that contains all the molecules within a structure.

Usage: module volume\_creation\_module

**MOVIE\_NAME** = ('movie\_name', <class 'str'>)

Name of the movie recorded by kwave.

Usage: adapter KwaveAcousticForwardModel

**NOISE\_FREQUENCY** = 'noise\_frequency'

Frequency of the noise model.

Usage: module core.processing\_components.noise

**NOISE\_MAX** = 'noise\_max'

Max of a noise model.

Usage: module core.processing\_components.noise

**NOISE\_MEAN** = 'noise\_mean'

Mean of a noise model.

Usage: module core.processing\_components.noise

**NOISE\_MIN** = 'noise\_min'

Min of a noise model.

Usage: module core.processing\_components.noise

**NOISE\_MODE** = 'noise\_mode'

The mode tag of a noise model is used to differentiate between

Tags.NOISE\_MODE\_ADDITIVE and Tags.NOISE\_MODE\_MULTIPLICATIVE.

Usage: module core.processing\_components.noise

**NOISE\_MODE\_ADDITIVE** = 'noise\_mode\_additive'

A noise model shall be applied additively  $s_n = s + n$ .

Usage: module core.processing\_components.noise

**NOISE\_MODE\_MULTIPLICATIVE** = 'noise\_mode\_multiplicative'

A noise model shall be applied multiplicatively  $s_n = s * n$ .

Usage: module core.processing\_components.noise

**NOISE\_NON\_NEGATIVITY\_CONSTRAINT** = 'noise\_non\_negativity\_constraint'

Defines if after the noise model negative values shall be allowed.

Usage: module core.processing\_components.noise

**NOISE\_SCALE** = 'noise\_scale'

Scale of a noise model.

Usage: module core.processing\_components.noise

**NOISE\_SHAPE** = 'noise\_shape'

Shape of a noise model.

Usage: module core.processing\_components.noise

**NOISE\_STD** = *'noise\_std'*

Standard deviation of a noise model.

Usage: module core.processing\_components.noise

**OPTICAL\_MODEL** = (*'optical\_model'*, <class 'str'>)

Choice of the used optical model.

Usage: module optical\_simulation\_module

**OPTICAL\_MODEL\_BINARY\_PATH** = (*'optical\_model\_binary\_path'*, <class 'str'>)

Absolute path of the location of the optical forward model binary.

Usage: module optical\_simulation\_module

**OPTICAL\_MODEL\_FLUENCE** = *'fluence'*

Name of the optical forward model output fluence field in the SIMPA output file.

Usage: naming convention

**OPTICAL\_MODEL\_ILLUMINATION\_GEOMETRY\_JSON\_FILE** = (*'optical\_model\_illumination\_geometry\_json\_file'*, <class 'str'>)

Absolute path of the location of the JSON file containing the IPASC-formatted optical forward model illumination geometry.

Usage: module optical\_simulation\_module

**OPTICAL\_MODEL\_INITIAL\_PRESSURE** = *'initial\_pressure'*

Name of the optical forward model output initial pressure field in the SIMPA output file.

Usage: naming convention

**OPTICAL\_MODEL\_MCX** = *'mcx'*

Corresponds to the mcx simulation.

Usage: module optical\_simulation\_module, naming convention

**OPTICAL\_MODEL\_NUMBER\_PHOTONS** = (*'optical\_model\_number\_of\_photons'*, (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Number of photons used in the optical simulation.

Usage: module optical\_simulation\_module

**OPTICAL\_MODEL\_OUTPUT\_NAME** = *'optical\_forward\_model\_output'*

Name of the optical forward model output field in the SIMPA output file.

Usage: naming convention

**OPTICAL\_MODEL\_SETTINGS** = (*'optical\_model\_settings'*, <class 'dict'>)

Optical model settings

**OPTICAL\_MODEL\_TEST** = *'simpa\_tests'*

Corresponds to an adapter for testing purposes only.

Usage: module optical\_simulation\_module, naming convention

**OPTICAL\_MODEL\_UNITS** = *'units'*

Name of the optical forward model output units field in the SIMPA output file.

Usage: naming convention

**ORIGINAL\_DATA** = *'original\_data'*

Name of the simulation outputs as original data in the SIMPA output file.

Usage: naming convention

**PARALLELEPIPED\_STRUCTURE** = *'ParallelepipedStructure'*

Corresponds to the ParallelepipedStructure in the structure\_library.

Usage: module volume\_creation\_module, naming\_convention

**PMLAlpha** = (*'pml\_alpha'*, (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Alpha coefficient of the “perfectly matched layer” (PML) around the simulated volume in kwave.

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter, naming convention

**PMLInside** = ('pml\_inside', (<class 'bool'>, <class 'bool'>, <class 'numpy.bool\_'>))

If True, the “perfectly matched layer” (PML) in kwave is located inside the volume.

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter, naming convention

**PMLSize** = ('pml\_size', (<class 'list'>, <class 'tuple'>, <class 'numpy.ndarray'>))

Size of the “perfectly matched layer” (PML) around the simulated volume in kwave.

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter, naming convention

**PRIORITY** = ('priority', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Number that corresponds to a priority of the assigned structure. If another structure occupies the same voxel in a volume, the structure with a higher priority will be preferred.

Usage: adapter versatile\_volume\_creator

**PROPERTY\_ABSORPTION\_PER\_CM** = 'mua'

Optical absorption of the generated volume/structure in 1/cm.

Usage: SIMPA package, naming convention

**PROPERTY\_ALPHA\_COEFF** = 'alpha\_coeff'

Acoustic attenuation of kwave of the generated volume/structure in dB/cm/MHz.

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter, naming convention

**PROPERTY\_ALPHA\_POWER** = ('medium\_alpha\_power', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Exponent of the exponential acoustic attenuation law of kwave.

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter, naming convention

**PROPERTY\_ANISOTROPY** = 'g'

Optical scattering anisotropy of the generated volume/structure.

Usage: SIMPA package, naming convention

**PROPERTY\_DENSITY** = 'density'

Density of the generated volume/structure in kg/m<sup>3</sup>.

Usage: SIMPA package, naming convention

**PROPERTY\_DIRECTIVITY\_ANGLE** = 'directivity\_angle'

Directionality of the sensors in kwave of the used PA device.

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter, naming convention

**PROPERTY\_GRUNEISEN\_PARAMETER** = 'gamma'

We define PROPERTY\_GRUNEISEN\_PARAMETER to contain all wavelength-independent constituents of the PA signal. This means that it contains the percentage of absorbed light converted into heat. Naturally, one could make an argument that this should not be the case, however, it simplifies the usage of this tool.

Usage: SIMPA package, naming convention

**PROPERTY\_INTRINSIC\_EULER\_ANGLE** = 'intrinsic\_euler\_angle'

Intrinsic euler angles of the detector elements in the kWaveArray.

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter, naming convention

**PROPERTY\_OXYGENATION** = 'oxy'

Oxygenation of the generated volume/structure.

Usage: SIMPA package, naming convention

**PROPERTY\_SCATTERING\_PER\_CM** = 'mus'

Optical scattering (NOT REDUCED SCATTERING mus'! mus'=mus\*(1-g) ) of the generated volume/structure in 1/cm.

Usage: SIMPA package, naming convention

**PROPERTY\_SEGMENTATION** = 'seg'

Segmentation of the generated volume/structure.

Usage: SIMPA package, naming convention

**PROPERTY\_SENSOR\_MASK** = *'sensor\_mask'*

Sensor mask of kwave of the used PA device.

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter, naming convention

**PROPERTY\_SPEED\_OF\_SOUND** = *'sos'*

Speed of sound of the generated volume/structure in m/s.

Usage: SIMPA package, naming convention

**PlotPML** = (*'plot\_pml'*, (*<class 'bool'>*, *<class 'bool'>*, *<class 'numpy.bool\_'>*))

If True, the “perfectly matched layer” (PML) around the simulated volume in kwave is plotted.

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter, naming convention

**RANDOM\_SEED** = (*'random\_seed'*, (*<class 'int'>*, *<class 'numpy.integer'>*))

Random seed for numpy and torch.

Usage: SIMPA package

**RECONSTRUCTED\_DATA** = *'reconstructed\_data'*

Name of the reconstructed data field in the SIMPA output file.

Usage: naming convention

**RECONSTRUCTED\_DATA\_NOISE** = *'reconstructed\_data\_noise'*

Name of the reconstructed data with applied noise field in the SIMPA output file.

Usage: naming convention

**RECONSTRUCTION\_ALGORITHM** = (*'reconstruction\_algorithm'*, *<class 'str'>*)

Choice of the used reconstruction algorithm.

Usage: module reconstruction\_module

**RECONSTRUCTION\_ALGORITHM\_DAS** = *'DAS'*

Corresponds to the reconstruction algorithm DAS with the MitkBeamformingAdapter.

Usage: module reconstruction\_module, naming convention

**RECONSTRUCTION\_ALGORITHM\_DMAS** = *'DMAS'*

Corresponds to the reconstruction algorithm DMAS with the MitkBeamformingAdapter.

Usage: module reconstruction\_module, naming convention

**RECONSTRUCTION\_ALGORITHM\_PYTORCH\_DAS** = *'PyTorch\_DAS'*

Corresponds to the reconstruction algorithm DAS with the PyTorchDASAdapter.

Usage: module reconstruction\_module, naming convention

**RECONSTRUCTION\_ALGORITHM\_SDMAS** = *'sDMAS'*

Corresponds to the reconstruction algorithm sDMAS with the MitkBeamformingAdapter.

Usage: module reconstruction\_module, naming convention

**RECONSTRUCTION\_ALGORITHM\_TEST** = *'TEST'*

Corresponds to an adapter for testing purposes only.

Usage: module reconstruction\_module, naming convention

**RECONSTRUCTION\_ALGORITHM\_TIME\_REVERSAL** = *'time\_reversal'*

Corresponds to the reconstruction algorithm Time Reversal with TimeReversalAdapter.

Usage: module reconstruction\_module, naming convention

**RECONSTRUCTION\_APODIZATION\_BOX** = *'BoxApodization'*

Corresponds to the box window function for apodization.

Usage: adapter PyTorchDASAdapter, naming convention

**RECONSTRUCTION\_APODIZATION\_HAMMING** = *'HammingApodization'*

Corresponds to the Hamming window function for apodization.

Usage: adapter PyTorchDASAdapter, naming convention

**RECONSTRUCTION\_APODIZATION\_HANN** = *'HannApodization'*

Corresponds to the Hann window function for apodization.

Usage: adapter PyTorchDASAdapter, naming convention

**RECONSTRUCTION\_APODIZATION\_METHOD** = ('reconstruction\_apodization\_method', <class 'str'>)

Choice of the apodization method used, i.e. window functions .

Usage: adapter PyTorchDASAdapter

**RECONSTRUCTION\_BMODE\_AFTER\_RECONSTRUCTION** = 'Envelope\_Detection\_after\_Reconstruction'

Specifies whether an envelope detection should be performed after reconstruction, default is False Usage:

adapter PyTorchDASAdapter, naming convention

**RECONSTRUCTION\_BMODE\_BEFORE\_RECONSTRUCTION** = 'Envelope\_Detection\_before\_Reconstruction'

Specifies whether an envelope detection should be performed before reconstruction, default is False Usage:

adapter PyTorchDASAdapter, naming convention

**RECONSTRUCTION\_BMODE\_METHOD** = ('reconstruction\_bmode\_method', <class 'str'>)

Choice of the B-Mode method used in the Mitk Beamforming.

Usage: adapter MitkBeamformingAdapter

**RECONSTRUCTION\_BMODE\_METHOD\_ABS** = 'Abs'

Corresponds to the absolute value as the B-Mode method used in the Mitk Beamforming.

Usage: adapter MitkBeamformingAdapter, naming convention

**RECONSTRUCTION\_BMODE\_METHOD\_HILBERT\_TRANSFORM** = 'EnvelopeDetection'

Corresponds to the Hilbert transform as the B-Mode method used in the Mitk Beamforming.

Usage: adapter MitkBeamformingAdapter, naming convention

**RECONSTRUCTION\_INVERSE\_CRIME** = ('reconstruction\_inverse\_crime', (<class 'bool'>, <class 'bool'>, <class 'numpy.bool\_'>))

If True, the Time Reversal reconstruction will commit the “inverse crime”.

Usage: TimeReversalAdapter

**RECONSTRUCTION\_MITK\_BINARY\_PATH** = ('reconstruction\_mitk\_binary\_path', <class 'str'>)

Absolute path to the Mitk Beamforming script.

Usage: adapter MitkBeamformingAdapter

**RECONSTRUCTION\_MITK\_SETTINGS\_XML** = ('reconstruction\_mitk\_settings\_xml', <class 'str'>)

Absolute path to the Mitk Beamforming script settings.

Usage: adapter MitkBeamformingAdapter

**RECONSTRUCTION\_MODE** = ('reconstruction\_mode', <class 'str'>)

Choice of the reconstruction mode used in the Backprojection.

Usage: adapter BackprojectionAdapter

**RECONSTRUCTION\_MODEL\_SETTINGS** = ('reconstruction\_model\_settings', <class 'dict'>)

” Reconstruction Model Settings

**RECONSTRUCTION\_MODE\_DIFFERENTIAL** = 'differential'

Corresponds to the differential mode used in the Backprojection.

Usage: adapter BackprojectionAdapter, naming\_convention

**RECONSTRUCTION\_MODE\_FULL** = 'full'

Corresponds to the full mode used in the Backprojection.

Usage: adapter BackprojectionAdapter, naming\_convention

**RECONSTRUCTION\_MODE\_PRESSURE** = 'pressure'

Corresponds to the pressure mode used in the Backprojection.

Usage: adapter BackprojectionAdapter, naming\_convention

**RECONSTRUCTION\_OUTPUT\_NAME** = ('reconstruction\_result', <class 'str'>)

Absolute path of the image reconstruction result.

Usage: adapter MitkBeamformingAdapter

**RECONSTRUCTION\_PERFORM\_BANDPASS\_FILTERING** = ('reconstruction\_perform\_bandpass\_filtering', (<class 'bool'>, <class 'bool'>, <class 'numpy.bool\_'>))

Whether bandpass filtering should be applied or not. Default should be True

Usage: adapter PyTorchDASAdapter

**RECORDMOVIE** = ('record\_movie', (<class 'bool'>, <class 'bool'>, <class 'numpy.bool\_'>))

If True, a movie of the kwave simulation will be recorded.

Usage: adapter KwaveAcousticForwardModel

**RECTANGULAR\_CUBOID\_STRUCTURE** = 'RectangularCuboidStructure'

Corresponds to the RectangularCuboidStructure in the structure\_library.

Usage: module volume\_creation\_module, naming\_convention

**SEGMENTATION\_CLASS\_MAPPING** = ('segmentation\_class\_mapping', <class 'dict'>)

Mapping that assigns every class in the INPUT\_SEGMENTATION\_VOLUME a MOLECULE\_COMPOSITION.

Usage: adapter segmentation\_based\_volume\_creator

**SENSOR\_BANDWIDTH\_PERCENT** = ('sensor\_bandwidth', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Sensor bandwidth in kwave.

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter, naming convention

**SENSOR\_CENTER\_FREQUENCY\_HZ** = ('sensor\_center\_frequency', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Sensor center frequency in kwave.

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter, naming convention

**SENSOR\_CONCAVE** = 'concave'

Indicates that the geometry of the used PA device in the Mitk Beamforming is concave.

Usage: adapter MitkBeamformingAdapter, naming convention

**SENSOR\_DIRECTIVITY\_PATTERN** = 'sensor\_directivity\_pattern'

Sensor directivity pattern of the sensor in kwave. Default should be "pressure".

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter, naming convention

**SENSOR\_DIRECTIVITY\_SIZE\_M** = ('sensor\_directivity\_size', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Size of each detector element in kwave.

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter, naming convention

**SENSOR\_ELEMENT\_POSITIONS** = 'sensor\_element\_positions'

Number of detector elements that fit into the generated volume if the dimensions and/or spacing of the generated volume were not highly resolved enough to be sufficient for the selected PA device.

Usage: module acoustic\_forward\_module, naming convention

**SENSOR\_LINEAR** = 'linear'

Indicates that the geometry of the used PA device in the Mitk Beamforming is linear.

Usage: adapter MitkBeamformingAdapter, naming convention

**SENSOR\_NUM\_ELEMENTS** = ('sensor\_num\_elements', (<class 'int'>, <class 'numpy.integer'>))

Number of detector elements for kwave if no device was selected.

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter, naming convention

**SENSOR\_NUM\_USED\_ELEMENTS** = ('sensor\_num\_used\_elements', (<class 'int'>, <class 'numpy.integer'>))

Number of detector elements that fit into the generated volume if the dimensions and/or spacing of the generated volume were not highly resolved enough to be sufficient for the selected PA device.

Usage: module acoustic\_forward\_module, naming convention

**SENSOR\_PITCH\_MM** = 'sensor\_pitch\_mm'

Pitch of detector elements of the used PA device.

Usage: adapter AcousticForwardModelKWaveAdapter, naming convention

**SENSOR\_RADIUS\_MM** = 'sensor\_radius\_mm'

Radius of a concave geometry of the used PA device.

Usage: adapter AcousticForwardModelKWaveAdapter, naming convention

**SENSOR\_RECORD** = ('sensor\_record', <class 'str'>)

Sensor Record mode of the sensor in kwave. Default should be "p".

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter, naming convention

**SENSOR\_SAMPLING\_RATE\_MHZ** = ('sensor\_sampling\_rate\_mhz', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Sampling rate of the used PA device.

Usage: adapter KwaveAcousticForwardModel, adapter TimeReversalAdapter, naming convention

**SETTINGS** = 'settings'

Location of the simulation settings in the SIMPA output file.

Usage: naming convention

**SETTINGS\_JSON** = ('settings\_json', (<class 'bool'>, <class 'numpy.bool\_'>))

If True, the SIMPA settings are saved in a .json file.

Usage: SIMPA package

**SETTINGS\_JSON\_PATH** = ('settings\_json\_path', <class 'str'>)

Absolute path to a .json file if SETTINGS\_JSON is set to True. Usage: SIMPA package

**SIMPA\_OUTPUT\_NAME** = 'simpa\_output.hdf5'

Default filename of the SIMPA output if not specified otherwise.

Usage: SIMPA package, naming convention

**SIMPA\_OUTPUT\_PATH** = ('simpa\_output\_path', <class 'str'>)

Default path of the SIMPA output if not specified otherwise.

Usage: SIMPA package

**SIMULATE\_DEFORMED\_LAYERS** = ('simulate\_deformed\_layers', <class 'bool'>)

If True, the horizontal layers are deformed according to the DEFORMED\_LAYERS\_SETTINGS.

Usage: adapter versatile\_volume\_creation

**SIMULATIONS** = 'simulations'

Location of the simulation outputs in the SIMPA output file.

Usage: naming convention

**SIMULATION\_PATH** = ('simulation\_path', <class 'str'>)

Absolute path to the folder where the SIMPA output is saved.

Usage: SIMPA package

**SIMULATION\_PROPERTIES** = 'simulation\_properties'

Location of the simulation properties in the SIMPA output file.

Usage: naming convention

**SPACING\_MM** = ('voxel\_spacing\_mm', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Isotropic extent of one voxels in mm in the generated volume.

Usage: SIMPA package

**SPHERICAL\_STRUCTURE** = 'SphericalStructure'

Corresponds to the SphericalStructure in the structure\_library.

Usage: module volume\_creation\_module, naming\_convention

**STRUCTURES** = ('structures', <class 'dict'>)

Settings dictionary which contains all the structures that should be generated inside the volume.

Usage: module volume\_creation\_module

**STRUCTURE\_BIFURCATION\_LENGTH\_MM** = ('structure\_bifurcation\_length\_mm', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Length after which a VesselStructure will bifurcate.

Usage: adapter versatile\_volume\_creation, class VesselStructure

**STRUCTURE\_CURVATURE\_FACTOR** = ('structure\_curvature\_factor', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Factor that determines how strongly a vessel tree is curved.

Usage: adapter versatile\_volume\_creation, class VesselStructure

**STRUCTURE\_DIRECTION** = ('structure\_direction', (<class 'list'>, <class 'tuple'>, <class 'numpy.ndarray'>))

Direction as [x, y, z] vector starting from STRUCTURE\_START\_MM in which the vessel will grow.

Usage: adapter versatile\_volume\_creation, class VesselStructure

**STRUCTURE\_ECCENTRICITY** = ('structure\_excentricity', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>, <class 'numpy.ndarray'>))

Eccentricity of the structure.

Usage: adapter versatile\_volume\_creation, class EllipticalTubularStructure

**STRUCTURE\_END\_MM** = ('structure\_end', (<class 'list'>, <class 'tuple'>, <class 'numpy.ndarray'>))

Ending of the structure as [x, y, z] coordinates in the generated volume.

Usage: adapter versatile\_volume\_creation, class GeometricalStructure

**STRUCTURE\_FIRST\_EDGE\_MM** = ('structure\_first\_edge\_mm', (<class 'list'>, <class 'tuple'>, <class 'numpy.ndarray'>))

Edge of the structure as [x, y, z] vector starting from STRUCTURE\_START\_MM in the generated volume.

Usage: adapter versatile\_volume\_creation, class ParallelepipedStructure

**STRUCTURE\_RADIUS\_MM** = ('structure\_radius', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>, <class 'numpy.ndarray'>))

Radius of the structure.

Usage: adapter versatile\_volume\_creation, class GeometricalStructure

**STRUCTURE\_RADIUS\_VARIATION\_FACTOR** = ('structure\_radius\_variation\_factor', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))

Factor that determines how strongly a the radius of vessel tree varies.

Usage: adapter versatile\_volume\_creation, class VesselStructure

**STRUCTURE\_SECOND\_EDGE\_MM** = ('structure\_second\_edge\_mm', (<class 'list'>, <class 'tuple'>, <class 'numpy.ndarray'>))

Edge of the structure as [x, y, z] vector starting from STRUCTURE\_START\_MM in the generated volume.

Usage: adapter versatile\_volume\_creation, class ParallelepipedStructure

**STRUCTURE\_SEGMENTATION\_TYPE** = 'structure\_segmentation\_type'

Defines the structure segmentation type to one segmentation type in SegmentationClasses.

Usage: module volume\_creation\_module, naming convention

**STRUCTURE\_START\_MM** = ('structure\_start', (<class 'list'>, <class 'tuple'>, <class 'numpy.ndarray'>))

Beginning of the structure as [x, y, z] coordinates in the generated volume.

Usage: adapter versatile\_volume\_creation, class GeometricalStructure

**STRUCTURE\_THIRD\_EDGE\_MM** = ('structure\_third\_edge\_mm', (<class 'list'>, <class 'tuple'>, <class 'numpy.ndarray'>))

Edge of the structure as [x, y, z] vector starting from STRUCTURE\_START\_MM in the generated volume.

Usage: adapter versatile\_volume\_creation, class ParallelepipedStructure

**STRUCTURE\_TYPE** = ('structure\_type', <class 'str'>)

Defines the structure type to one structure in the structure\_library.

Usage: module volume\_creation\_module



**STRUCTURE\_X\_EXTENT\_MM** = ('structure\_x\_extent\_mm', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))  
X-extent of the structure in the generated volume.  
Usage: adapter versatile\_volume\_creation, class RectangularCuboidStructure

**STRUCTURE\_Y\_EXTENT\_MM** = ('structure\_y\_extent\_mm', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))  
Y-extent of the structure in the generated volume.  
Usage: adapter versatile\_volume\_creation, class RectangularCuboidStructure

**STRUCTURE\_Z\_EXTENT\_MM** = ('structure\_z\_extent\_mm', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))  
Z-extent of the structure in the generated volume.  
Usage: adapter versatile\_volume\_creation, class RectangularCuboidStructure

**TIME\_SERIES\_DATA** = 'time\_series\_data'  
Name of the time series data field in the SIMPA output file.  
Usage: naming convention

**TIME\_SERIES\_DATA\_NOISE** = 'time\_series\_data\_noise'  
Name of the time series data with applied noise field in the SIMPA output file.  
Usage: naming convention

**TIME\_STEP** = ('time\_step', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))  
Temporal resolution of mcx.  
Usage: adapter mcx\_adapter

**TISSUE\_PROPERTIES\_OUPUT\_NAME** = 'properties'  
Name of the simulation properties field in the SIMPA output file.  
Usage: naming convention

**TOTAL\_TIME** = ('total\_time', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))  
Total simulated time in mcx.  
Usage: adapter mcx\_adapter

**TUKEY\_WINDOW\_ALPHA** = ('tukey\_window\_alpha', (<class 'int'>, <class 'numpy.integer'>, <class 'float'>, <class 'float'>))  
Sets alpha value of Tukey window between 0 (similar to box window) and 1 (similar to Hann window). Default is 0.5  
Usage: adapter PyTorchDASAdapter

**UNITS\_ARBITRARY** = 'arbitrary\_unity'  
Define arbitrary units if no units were given in the settings.  
Usage: module optical\_simulation\_module, naming convention

**UNITS\_PRESSURE** = 'newton\_per\_meters\_squared'  
Standard units used in the SIMPA framework.  
Usage: module optical\_simulation\_module, naming convention

**UPSAMPLED\_DATA** = 'upsampled\_data'  
Name of the simulation outputs as upsampled data in the SIMPA output file.  
Usage: naming convention

**US\_GEL** = ('us\_gel', <class 'bool'>)  
If True, us gel is placed between the PA device and the simulated volume.  
Usage: SIMPA package

**VOLUME\_CREATION\_MODEL\_SETTINGS** = ('volume\_creation\_model\_settings', <class 'dict'>)  
" Volume Creation Model Settings

**VOLUME\_CREATOR** = ('volume\_creator', <class 'str'>)

Choice of the volume creator adapter.

Usage: module volume\_creation\_module, module device\_digital\_twins

**VOLUME\_CREATOR\_SEGMENTATION\_BASED** = 'volume\_creator\_segmentation\_based'

Corresponds to the SegmentationBasedVolumeCreator.

Usage: module volume\_creation\_module, naming convention

**VOLUME\_CREATOR\_VERSATILE** = 'volume\_creator\_versatile'

Corresponds to the ModelBasedVolumeCreator.

Usage: module volume\_creation\_module, naming convention

**VOLUME\_NAME** = ('volume\_name', <class 'str'>)

Name of the SIMPA output file.

Usage: SIMPA package

**WAVELENGTH** = ('wavelength', (<class 'int'>, <class 'numpy.integer'>))

Single wavelength used for the current simulation.

Usage: SIMPA package

**WAVELENGTHS** = ('wavelengths', (<class 'list'>, <class 'range'>, <class 'tuple'>, <class 'numpy.ndarray'>))

Iterable of all the wavelengths used for the simulation.

Usage: SIMPA package

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI)  
SPDX-License-Identifier: MIT

`class simpa.utils.tissue_properties.TissueProperties`

## Libraries

Another important aspect of the utils class is the libraries that are being provided. These contain compilations of literature values for the acoustic and optical properties of commonly used tissue.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI)  
SPDX-License-Identifier: MIT

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI)  
SPDX-License-Identifier: MIT

`class` `simpa.utils.libraries.molecule_library.MolecularComposition`  
(segmentation\_type=None, molecular\_composition\_settings=None)

**update\_internal\_properties ()**  
FIXME

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI)  
SPDX-License-Identifier: MIT

`class simpa.utils.libraries.literature_values.MorphologicalTissueProperties`

This class contains a listing of morphological tissue parameters as reported in literature. The listing is not the result of a meta analysis, but rather uses the best fitting paper at the time of implementation. Each of the fields is annotated with a literature reference or a descriptions of how the particular values were derived for tissue modelling.

`class simpa.utils.libraries.literature_values.OpticalTissueProperties`

This class contains a listing of optical tissue parameters as reported in literature. The listing is not the result of a meta analysis, but rather uses the best fitting paper at the time of implementation. Each of the fields is annotated with a literature reference or a descriptions of how the particular values were derived for tissue modelling.

`class simpa.utils.libraries.literature_values.StandardProperties`

This class contains a listing of default parameters that can be used. These values are sensible default values but are generally not backed up by proper scientific references, or are rather specific for internal use cases.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.utils.libraries.spectra_library.Spectrum (spectrum_name: str, wavelengths:
numpy.ndarray, values: numpy.ndarray)
```

An instance of this class represents the absorption spectrum over wavelength for a particular

```
get_value_for_wavelength (wavelength: int) → float
```

**Parameters:** **wavelength** – the wavelength to retrieve a optical absorption value for [ $\text{cm}^{-1}$ ]. Must be an integer value between the minimum and maximum wavelength.

**Returns:** the best matching linearly interpolated absorption value for the given wavelength.

```
get_value_over_wavelength ()
```

**Returns:** numpy array with the available wavelengths and the corresponding absorption properties

```
simpa.utils.libraries.spectra_library.view_absorption_spectra (save_path=None)
```

Opens a matplotlib plot and visualizes the available absorption spectra.

**Parameters:** **save\_path** – If not None, then the figure will be saved as a png file to the destination.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.utils.libraries.tissue_library.MolecularCompositionGenerator
```

The MolecularCompositionGenerator is a helper class to facilitate the creation of a MolecularComposition instance.

```
class simpa.utils.libraries.tissue_library.TissueLibrary
TODO
```

```
blood (oxygenation=None)
```

**Returns:** a settings dictionary containing all min and max parameters fitting for full blood.

```
bone ()
```

**Returns:** a settings dictionary containing all min and max parameters fitting for full blood.

```
constant (mua=1e-10, mus=1e-10, g=1e-10)
TODO
```

```
dermis (background_oxy=None, blood_volume_fraction=None)
```

**Returns:** a settings dictionary containing all min and max parameters fitting for dermis tissue.

```
epidermis (melanosom_volume_fraction=None)
```

**Returns:** a settings dictionary containing all min and max parameters fitting for epidermis tissue.

```
get_blood_volume_fractions (total_blood_volume_fraction=1e-10, oxygenation=1e-10)
TODO
```

```
muscle (background_oxy=None, blood_volume_fraction=None)
```

**Returns:** a settings dictionary containing all min and max parameters fitting for generic background tissue.

```
subcutaneous_fat (background_oxy=0.5)
```

**Returns:** a settings dictionary containing all min and max parameters fitting for subcutaneous fat tissue.

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.utils.libraries.structure_library.Background (global_settings:
simpa.utils.settings.Settings, background_settings: simpa.utils.settings.Settings =
None)
```

Defines a background that fills the whole simulation volume. It is always given the priority of 0 so that other structures can overwrite it when necessary.

Example usage:

```
background_dictionary = Settings()
background_dictionary[Tags.MOLECULE_COMPOSITION] = TISSUE_LIBRARY.constant(0.1, 100.0, 0.9)
background_dictionary[Tags.STRUCTURE_TYPE] = Tags.BACKGROUND
```

**get\_enclosed\_indices ()**

Gets indices of the voxels that are either entirely or partially occupied by the GeometricalStructure. :return: mask for a numpy array

**get\_params\_from\_settings (single\_structure\_settings)**

Gets all the parameters required for the specific GeometricalStructure. :param single\_structure\_settings: Settings which describe the specific GeometricalStructure. :return: Tuple of parameters

**to\_settings ()** → dict

Creates a Settings dictionary which contains all the parameters needed to create the same GeometricalStructure again. :return: A tuple containing the settings key and the needed entries

```
class simpa.utils.libraries.structure_library.CircularTubularStructure (global_settings:
simpa.utils.settings.Settings, single_structure_settings:
simpa.utils.settings.Settings = None)
```

Defines a circular tube which is defined by a start and end point as well as a radius. This structure implements partial volume effects. The tube can be set to adhere to a deformation defined by the simpa.utils.deformation\_manager. The start and end points of the tube will then be shifted along the z-axis accordingly. Example usage:

```
# single_structure_settings initialization
structure = Settings()
structure[Tags.PRIORITY] = 9
structure[Tags.STRUCTURE_START_MM] = [50, 0, 50]
structure[Tags.STRUCTURE_END_MM] = [50, 100, 50]
structure[Tags.STRUCTURE_RADIUS_MM] = 5
structure[Tags.MOLECULE_COMPOSITION] = TISSUE_LIBRARY.blood()
structure[Tags.CONSIDER_PARTIAL_VOLUME] = True
structure[Tags.ADHERE_TO_DEFORMATION] = True
structure[Tags.STRUCTURE_TYPE] = Tags.CIRCULAR_TUBULAR_STRUCTURE
```

**get\_enclosed\_indices ()**

Gets indices of the voxels that are either entirely or partially occupied by the GeometricalStructure. :return: mask for a numpy array

**get\_params\_from\_settings (single\_structure\_settings)**

Gets all the parameters required for the specific GeometricalStructure. :param single\_structure\_settings: Settings which describe the specific GeometricalStructure. :return: Tuple of parameters

**to\_settings ()**

Creates a Settings dictionary which contains all the parameters needed to create the same GeometricalStructure again. :return: A tuple containing the settings key and the needed entries

```
class simpa.utils.libraries.structure_library.EllipticalTubularStructure
(global_settings: simpa.utils.settings.Settings, single_structure_settings:
simpa.utils.settings.Settings = None)
```

Defines an elliptical tube which is defined by a start and end point as well as a radius and an eccentricity. The elliptical geometry corresponds to a circular tube of the specified radius which is compressed along the z-axis until it reaches the specified eccentricity under the assumption of a constant volume. This structure implements partial

volume effects. The tube can be set to adhere to a deformation defined by the `simpa.utils.deformation_manager`. The start and end points of the tube will then be shifted along the z-axis accordingly. Example usage:

```
# single_structure_settings initialization
structure = Settings()

structure[Tags.PRIORITY] = 9
structure[Tags.STRUCTURE_START_MM] = [50, 0, 50]
structure[Tags.STRUCTURE_END_MM] = [50, 100, 50]
structure[Tags.STRUCTURE_RADIUS_MM] = 5
structure[Tags.STRUCTURE_ECCENTRICITY] = 0.8
structure[Tags.MOLECULE_COMPOSITION] = TISSUE_LIBRARY.blood()
structure[Tags.CONSIDER_PARTIAL_VOLUME] = True
structure[Tags.ADHERE_TO_DEFORMATION] = True
structure[Tags.STRUCTURE_TYPE] = Tags.ELLIPTICAL_TUBULAR_STRUCTURE
```

#### **get\_enclosed\_indices ()**

Gets indices of the voxels that are either entirely or partially occupied by the `GeometricalStructure`. :return: mask for a numpy array

#### **get\_params\_from\_settings (single\_structure\_settings)**

Gets all the parameters required for the specific `GeometricalStructure`. :param single\_structure\_settings: Settings which describe the specific `GeometricalStructure`. :return: Tuple of parameters

#### **to\_settings ()**

Creates a `Settings` dictionary which contains all the parameters needed to create the same `GeometricalStructure` again. :return : A tuple containing the settings key and the needed entries

```
class simpa.utils.libraries.structure_library.GeometricalStructure (global_settings:
simpa.utils.settings.Settings, single_structure_settings:
simpa.utils.settings.Settings = None)
```

Base class for all model-based structures for `ModelBasedVolumeCreator`. A `GeometricalStructure` has an internal representation of its own geometry. This is represented by `self.geometrical_volume` which is a 3D array that defines for every voxel within the simulation volume if it is enclosed in the `GeometricalStructure` or if it is outside. Most of the `GeometricalStructures` implement a partial volume effect. So if a voxel has the value 1, it is completely enclosed by the `GeometricalStructure`. If a voxel has a value between 0 and 1, that fraction of the volume is occupied by the `GeometricalStructure`. If a voxel has the value 0, it is outside of the `GeometricalStructure`.

#### **fill\_internal\_volume ()**

Fills `self.geometrical_volume` of the `GeometricalStructure`.

#### **abstract get\_enclosed\_indices ()**

Gets indices of the voxels that are either entirely or partially occupied by the `GeometricalStructure`. :return: mask for a numpy array

#### **abstract get\_params\_from\_settings (single\_structure\_settings)**

Gets all the parameters required for the specific `GeometricalStructure`. :param single\_structure\_settings: Settings which describe the specific `GeometricalStructure`. :return: Tuple of parameters

#### **properties\_for\_wavelength (wavelength) → simpa.utils.tissue\_properties.TissueProperties**

Returns the values corresponding to each optical/acoustic property used in SIMPA. :param wavelength: Wavelength of the queried properties :return: optical/acoustic properties

#### **abstract to\_settings () → simpa.utils.settings.Settings**

Creates a `Settings` dictionary which contains all the parameters needed to create the same `GeometricalStructure` again. :return : A tuple containing the settings key and the needed entries

```
class simpa.utils.libraries.structure_library.HorizontalLayerStructure (global_settings:
simpa.utils.settings.Settings, single_structure_settings:
simpa.utils.settings.Settings = None)
```

Defines a `Layer` structure which spans the xy-plane in the SIMPA axis convention. The thickness of the layer is defined along the z-axis. This layer can be deformed by the `simpa.utils.deformation_manager`. Example usage:

```
# single_structure_settings initialization
structure = Settings()

structure[Tags.PRIORITY] = 10
structure[Tags.STRUCTURE_START_MM] = [0, 0, 0]
structure[Tags.STRUCTURE_END_MM] = [0, 0, 100]
structure[Tags.MOLECULE_COMPOSITION] =
```

```

TISSUE_LIBRARY.epidermis()      structure[Tags.CONSIDER_PARTIAL_VOLUME] = True
structure[Tags.ADHERE_TO_DEFORMATION] = True structure[Tags.STRUCTURE_TYPE] =
Tags.HORIZONTAL_LAYER_STRUCTURE

```

#### **get\_enclosed\_indices ()**

Gets indices of the voxels that are either entirely or partially occupied by the GeometricalStructure. :return: mask for a numpy array

#### **get\_params\_from\_settings (single\_structure\_settings)**

Gets all the parameters required for the specific GeometricalStructure. :param single\_structure\_settings: Settings which describe the specific GeometricalStructure. :return: Tuple of parameters

#### **to\_settings ()**

Creates a Settings dictionary which contains all the parameters needed to create the same GeometricalStructure again. :return : A tuple containing the settings key and the needed entries

```

class simpa.utils.libraries.structure_library.ParallelepipedStructure (global_settings:
simpa.utils.settings.Settings, single_structure_settings:
simpa.utils.settings.Settings = None)

```

Defines a parallelepiped which is defined by a start point and three edge vectors which originate from the start point. This structure currently does not implement partial volume effects. Example usage:

```

# single_structure_settings initialization structure = Settings()

structure[Tags.PRIORITY] = 9 structure[Tags.STRUCTURE_START_MM] = [25, 25, 25]
structure[Tags.STRUCTURE_FIRST_EDGE_MM] = [5, 1, 1]
structure[Tags.STRUCTURE_SECOND_EDGE_MM] = [1, 5, 1]
structure[Tags.STRUCTURE_THIRD_EDGE_MM] = [1, 1, 5] structure[Tags.MOLECULE_COMPOSITION] =
TISSUE_LIBRARY.muscle() structure[Tags.STRUCTURE_TYPE] = Tags.PARALLELEPIPED_STRUCTURE

```

#### **get\_enclosed\_indices ()**

Gets indices of the voxels that are either entirely or partially occupied by the GeometricalStructure. :return: mask for a numpy array

#### **get\_params\_from\_settings (single\_structure\_settings)**

Gets all the parameters required for the specific GeometricalStructure. :param single\_structure\_settings: Settings which describe the specific GeometricalStructure. :return: Tuple of parameters

#### **to\_settings ()**

Creates a Settings dictionary which contains all the parameters needed to create the same GeometricalStructure again. :return : A tuple containing the settings key and the needed entries

```

class simpa.utils.libraries.structure_library.RectangularCuboidStructure
(global_settings: simpa.utils.settings.Settings, single_structure_settings:
simpa.utils.settings.Settings = None)

```

Defines a rectangular cuboid (box) which is defined by a start point its extent along the x-, y-, and z-axis. This structure implements partial volume effects. The box can be set to adhere to a deformation defined by the simpa.utils.deformation\_manager. The start point of the box will then be shifted along the z-axis accordingly. Example usage:

```

# single_structure_settings initialization structure = Settings()

structure[Tags.PRIORITY] = 9 structure[Tags.STRUCTURE_START_MM] = [25, 25, 25]
structure[Tags.STRUCTURE_X_EXTENT_MM] = 40 structure[Tags.STRUCTURE_Y_EXTENT_MM] = 50
structure[Tags.STRUCTURE_Z_EXTENT_MM] = 60 structure[Tags.MOLECULE_COMPOSITION] =
TISSUE_LIBRARY.muscle() structure[Tags.CONSIDER_PARTIAL_VOLUME] = True
structure[Tags.ADHERE_TO_DEFORMATION] = True structure[Tags.STRUCTURE_TYPE] =
Tags.RECTANGULAR_CUBOID_STRUCTURE

```

#### **get\_enclosed\_indices ()**

Gets indices of the voxels that are either entirely or partially occupied by the GeometricalStructure. :return: mask for a numpy array

**get\_params\_from\_settings (single\_structure\_settings)**

Gets all the parameters required for the specific GeometricalStructure. :param single\_structure\_settings: Settings which describe the specific GeometricalStructure. :return: Tuple of parameters

**to\_settings ()**

Creates a Settings dictionary which contains all the parameters needed to create the same GeometricalStructure again. :return : A tuple containing the settings key and the needed entries

```
class      simpa.utils.libraries.structure_library.SphericalStructure      (global_settings:
simpa.utils.settings.Settings,                                             single_structure_settings:
simpa.utils.settings.Settings = None)
```

Defines a sphere which is defined by a start point and a radius. This structure implements partial volume effects. The sphere can be set to adhere to a deformation defined by the simpa.utils.deformation\_manager. The start point of the sphere will then be shifted along the z-axis accordingly. Example usage:

```
# single_structure_settings initialization structure = Settings()
```

```
structure[Tags.PRIORITY]      = 9      structure[Tags.STRUCTURE_START_MM] = [50, 50, 50]
structure[Tags.STRUCTURE_RADIUS_MM] = 10 structure[Tags.MOLECULE_COMPOSITION] =
TISSUE_LIBRARY.blood()      structure[Tags.CONSIDER_PARTIAL_VOLUME]      = True
structure[Tags.ADHERE_TO_DEFORMATION] = True structure[Tags.STRUCTURE_TYPE] =
Tags.SPHERICAL_STRUCTURE
```

**get\_enclosed\_indices ()**

Gets indices of the voxels that are either entirely or partially occupied by the GeometricalStructure. :return: mask for a numpy array

**get\_params\_from\_settings (single\_structure\_settings)**

Gets all the parameters required for the specific GeometricalStructure. :param single\_structure\_settings: Settings which describe the specific GeometricalStructure. :return: Tuple of parameters

**to\_settings ()**

Creates a Settings dictionary which contains all the parameters needed to create the same GeometricalStructure again. :return : A tuple containing the settings key and the needed entries

```
class      simpa.utils.libraries.structure_library.Structures      (settings:
simpa.utils.settings.Settings, volume_creator_settings: dict)
TODO
```

```
class      simpa.utils.libraries.structure_library.VesselStructure      (global_settings:
simpa.utils.settings.Settings,                                             single_structure_settings:
simpa.utils.settings.Settings = None)
```

Defines a vessel tree that is generated randomly in the simulation volume. The generation process begins at the start with a specified radius. The vessel grows roughly in the specified direction. The deviation is specified by the curvature factor. Furthermore, the radius of the vessel can vary depending on the specified radius variation factor. The bifurcation length defines how long a vessel can get until it will bifurcate. This structure implements partial volume effects. Example usage:

```
# single_structure_settings initialization structure_settings = Settings()
```

```
structure_settings[Tags.PRIORITY] = 10 structure_settings[Tags.STRUCTURE_START_MM] = [50, 0, 50]
structure_settings[Tags.STRUCTURE_DIRECTION] = [0, 1, 0]
structure_settings[Tags.STRUCTURE_RADIUS_MM] = 4
structure_settings[Tags.STRUCTURE_CURVATURE_FACTOR] = 0.05
structure_settings[Tags.STRUCTURE_RADIUS_VARIATION_FACTOR] = 1
structure_settings[Tags.STRUCTURE_BIFURCATION_LENGTH_MM] = 70
structure_settings[Tags.MOLECULE_COMPOSITION] = TISSUE_LIBRARY.blood()
structure_settings[Tags.CONSIDER_PARTIAL_VOLUME] = True
structure_settings[Tags.STRUCTURE_TYPE] = Tags.VESSEL_STRUCTURE
```

**fill\_internal\_volume ()**

Fills self.geometrical\_volume of the GeometricalStructure.

### **get\_enclosed\_indices ()**

Gets indices of the voxels that are either entirely or partially occupied by the GeometricalStructure. :return: mask for a numpy array

### **get\_params\_from\_settings (single\_structure\_settings)**

Gets all the parameters required for the specific GeometricalStructure. :param single\_structure\_settings: Settings which describe the specific GeometricalStructure. :return: Tuple of parameters

### **to\_settings ()**

Creates a Settings dictionary which contains all the parameters needed to create the same GeometricalStructure again. :return : A tuple containing the settings key and the needed entries

## Module: io\_handling

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
simpa.io_handling.io_hdf5.load_hdf5 (file_path, file_dictionary_path='/')
```

Loads a dictionary from an hdf5 file.

#### **Parameters:**

- **file\_path** – Path of the file to load the dictionary from.
- **file\_dictionary\_path** – Path in dictionary structure of hdf5 file to load the dictionary in.

**Returns:** Dictionary

```
simpa.io_handling.io_hdf5.save_hdf5 (save_item, file_path: str, file_dictionary_path: str = '/', file_compression: str = None)
```

Saves a dictionary with arbitrary content or an item of any kind to an hdf5-file with given filepath.

#### **Parameters:**

- **save\_item** – Dictionary to save.
- **file\_path** – Path of the file to save the dictionary in.
- **file\_dictionary\_path** – Path in dictionary structure of existing hdf5 file to store the dictionary in.
- **file\_compression** – possible file compression for the hdf5 output file. Values are: gzip, lzf and szip.

**Returns:** **Null**

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT

```
class simpa.io_handling.serialization.SIMPAJSONSerializer
```

TODO

```
def default (_object: object)
```

TODO

```
class simpa.io_handling.serialization.SIMPASerializer
```

TODO

```
def serialize (_object: object)
```

## Module: log

SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI) SPDX-License-Identifier: MIT



SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ  
 SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI)  
 SPDX-License-Identifier: MIT

`class simpa.log.file_logger.Logger`

The SIMPA Logger. The purpose of this class is to guarantee that the logging Config has been set and that logging strings are written to the same file throughout the entire simulation pipeline. Per default, the log file is located in the home directory as defined by `Path.home()`.

The log levels are defined the same way they are in the python logging module: **DEBUG**: Detailed information, typically of interest only when diagnosing problems. **INFO**: Confirmation that things are working as expected. **WARNING**: An indication that something unexpected happened, or indicative of some problem in the near future (e.g. 'disk space low'). The software is still working as expected. **ERROR**: Due to a more serious problem, the software has not been able to perform some function. **CRITICAL**: A serious error, indicating that the program itself may be unable to continue running.

**critical** (`msg`)

Logs a critical message to the logging system. **CRITICAL**: A serious error, indicating that the program itself may be unable to continue running.

**Parameters:** `msg` – the message to log

**debug** (`msg`)

Logs a debug message to the logging system. **DEBUG**: Detailed information, typically of interest only when diagnosing problems.

**Parameters:** `msg` – the message to log

**error** (`msg`)

Logs an error message to the logging system. **ERROR**: Due to a more serious problem, the software has not been able to perform some function.

**Parameters:** `msg` – the message to log

**info** (`msg`)

Logs an info message to the logging system. **INFO**: Confirmation that things are working as expected.

**Parameters:** `msg` – the message to log

**warning** (`msg`)

Logs a warning message to the logging system. **WARNING**: An indication that something unexpected happened, or indicative of some problem in the near future (e.g. 'disk space low'). The software is still working as expected.

**Parameters:** `msg` – the message to log

## Examples

### Performing an optical forward simulation

The file can be found in `simpa_examples/minimal_optical_simulation.py`:

```
"""
SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI)
SPDX-License-Identifier: MIT
"""

from simpa.utils import Tags, TISSUE_LIBRARY
from simpa.core.simulation import simulate
from simpa.utils.settings import Settings
from simpa.visualisation.matplotlib_data_visualisation import visualise_data
```

```

from simpa.core.device_digital_twins import *
import numpy as np
from simpa.core import VolumeCreationModelModelBasedAdapter, OpticalForwardModelMcxAdapter,
    GaussianNoiseProcessingComponent

from simpa.utils.path_manager import PathManager

# FIXME temporary workaround for newest Intel architectures
import os
os.environ["KMP_DUPLICATE_LIB_OK"] = "TRUE"

# TODO: Please make sure that a valid path_config.env file is located in your home directory
# point to the correct file in the PathManager().
path_manager = PathManager()

VOLUME_TRANSDUCER_DIM_IN_MM = 60
VOLUME_PLANAR_DIM_IN_MM = 30
VOLUME_HEIGHT_IN_MM = 60
SPACING = 0.5
RANDOM_SEED = 471
VOLUME_NAME = "MyVolumeName_" + str(RANDOM_SEED)

# If VISUALIZE is set to True, the simulation result will be plotted
VISUALIZE = True

def create_example_tissue():
    """
    This is a very simple example script of how to create a tissue definition.
    It contains a muscular background, an epidermis layer on top of the muscles
    and a blood vessel.
    """
    background_dictionary = Settings()
    background_dictionary[Tags.MOLECULE_COMPOSITION] = TISSUE_LIBRARY.constant(1e-4, 1e-4, 0)
    background_dictionary[Tags.STRUCTURE_TYPE] = Tags.BACKGROUND

    muscle_dictionary = Settings()
    muscle_dictionary[Tags.PRIORITY] = 1
    muscle_dictionary[Tags.STRUCTURE_START_MM] = [0, 0, 10]
    muscle_dictionary[Tags.STRUCTURE_END_MM] = [0, 0, 100]
    muscle_dictionary[Tags.MOLECULE_COMPOSITION] = TISSUE_LIBRARY.muscle()
    muscle_dictionary[Tags.CONSIDER_PARTIAL_VOLUME] = True
    muscle_dictionary[Tags.ADHERE_TO_DEFORMATION] = True
    muscle_dictionary[Tags.STRUCTURE_TYPE] = Tags.HORIZONTAL_LAYER_STRUCTURE

    vessel_1_dictionary = Settings()
    vessel_1_dictionary[Tags.PRIORITY] = 3
    vessel_1_dictionary[Tags.STRUCTURE_START_MM] = [VOLUME_TRANSDUCER_DIM_IN_MM/2,
                                                    10,
                                                    VOLUME_HEIGHT_IN_MM/2]
    vessel_1_dictionary[Tags.STRUCTURE_END_MM] = [VOLUME_TRANSDUCER_DIM_IN_MM/2,
                                                    12,
                                                    VOLUME_HEIGHT_IN_MM/2]
    vessel_1_dictionary[Tags.STRUCTURE_RADIUS_MM] = 3
    vessel_1_dictionary[Tags.MOLECULE_COMPOSITION] = TISSUE_LIBRARY.blood()
    vessel_1_dictionary[Tags.CONSIDER_PARTIAL_VOLUME] = True
    vessel_1_dictionary[Tags.STRUCTURE_TYPE] = Tags.CIRCULAR_TUBULAR_STRUCTURE

    epidermis_dictionary = Settings()
    epidermis_dictionary[Tags.PRIORITY] = 8

```

```

epidermis_dictionary[Tags.STRUCTURE_START_MM] = [0, 0, 9]
epidermis_dictionary[Tags.STRUCTURE_END_MM] = [0, 0, 10]
epidermis_dictionary[Tags.MOLECULE_COMPOSITION] = TISSUE_LIBRARY.epidermis()
epidermis_dictionary[Tags.CONSIDER_PARTIAL_VOLUME] = True
epidermis_dictionary[Tags.ADHERE_TO_DEFORMATION] = True
epidermis_dictionary[Tags.STRUCTURE_TYPE] = Tags.HORIZONTAL_LAYER_STRUCTURE

tissue_dict = Settings()
tissue_dict[Tags.BACKGROUND] = background_dictionary
tissue_dict["muscle"] = muscle_dictionary
tissue_dict["epidermis"] = epidermis_dictionary
tissue_dict["vessel_1"] = vessel_1_dictionary
return tissue_dict

# Seed the numpy random configuration prior to creating the global_settings file in
# order to ensure that the same volume
# is generated with the same random seed every time.

np.random.seed(RANDOM_SEED)

general_settings = {
    # These parameters set the general properties of the simulated volume
    Tags.RANDOM_SEED: RANDOM_SEED,
    Tags.VOLUME_NAME: VOLUME_NAME,
    Tags.SIMULATION_PATH: path_manager.get_hdf5_file_save_path(),
    Tags.SPACING_MM: SPACING,
    Tags.DIM_VOLUME_Z_MM: VOLUME_HEIGHT_IN_MM,
    Tags.DIM_VOLUME_X_MM: VOLUME_TRANSDUCER_DIM_IN_MM,
    Tags.DIM_VOLUME_Y_MM: VOLUME_PLANAR_DIM_IN_MM,
    Tags.WAVELENGTHS: [798],
    Tags.DIGITAL_DEVICE_POSITION: [VOLUME_TRANSDUCER_DIM_IN_MM/2,
                                   VOLUME_PLANAR_DIM_IN_MM/2,
                                   0],
    Tags.LOAD_AND_SAVE_HDF5_FILE_AT_THE_END_OF_SIMULATION_TO_MINIMISE_FILESIZE: True
}

settings = Settings(general_settings)

settings.set_volume_creation_settings({
    Tags.SIMULATE_DEFORMED_LAYERS: True,
    Tags.STRUCTURES: create_example_tissue()
})

settings.set_optical_settings({
    Tags.OPTICAL_MODEL_NUMBER_PHOTONS: 1e7,
    Tags.OPTICAL_MODEL_BINARY_PATH: path_manager.get_mcx_binary_path(),
    Tags.OPTICAL_MODEL: Tags.OPTICAL_MODEL_MCX,
    Tags.ILLUMINATION_TYPE: Tags.ILLUMINATION_TYPE_PENCIL,
    Tags.LASER_PULSE_ENERGY_IN_MILLIJOULE: 50
})

settings["noise_model_1"] = {
    Tags.NOISE_MEAN: 1.0,
    Tags.NOISE_STD: 0.1,
    Tags.NOISE_MODE: Tags.NOISE_MODE_MULTIPLICATIVE,
    Tags.DATA_FIELD: Tags.OPTICAL_MODEL_INITIAL_PRESSURE,
    Tags.NOISE_NON_NEGATIVITY_CONSTRAINT: True
}

pipeline = [
    VolumeCreationModelModelBasedAdapter(settings),
    OpticalForwardModelMcxAdapter(settings),

```

```

    GaussianNoiseProcessingComponent(settings, "noise_model_1")
]

class ExampleDeviceSlitIlluminationLinearDetector(PhotoacousticDevice):
    """
    This class represents a digital twin of a PA device with a slit as illumination next to
    """

    def __init__(self):
        super().__init__(device_position_mm=np.asarray([VOLUME_TRANSDUCER_DIM_IN_MM/2,
                                                         VOLUME_PLANAR_DIM_IN_MM/2, 0]))
        self.set_detection_geometry(LinearArrayDetectionGeometry())
        self.add_illumination_geometry(SlitIlluminationGeometry(slit_vector_mm=[20, 0, 0],
                                                                direction_vector_mm=[0, 0, 5]))

device = ExampleDeviceSlitIlluminationLinearDetector()

simulate(pipeline, settings, device)

if Tags.WAVELENGTH in settings:
    WAVELENGTH = settings[Tags.WAVELENGTH]
else:
    WAVELENGTH = 700

if VISUALIZE:
    visualise_data(path_manager.get_hdf5_file_save_path() + "/" + VOLUME_NAME + ".hdf5", WAVELENGTH,
                  log_scale=True)

```

## Performing a complete forward simulation with acoustic modeling, optical modeling, as well as image reconstruction

The file can be found in `simpa_examples/optical_and_acoustic_simulation.py`:

```

"""
SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI)
SPDX-License-Identifier: MIT
"""

from simpa.utils import Tags, TISSUE_LIBRARY

from simpa.core.simulation import simulate
from simpa.utils.settings import Settings
from simpa.visualisation.matplotlib_data_visualisation import visualise_data
import numpy as np
from simpa.utils.path_manager import PathManager
from simpa.core import ImageReconstructionModuleDelayAndSumAdapter, GaussianNoiseProcessingComponent,
    OpticalForwardModelMcxAdapter, AcousticForwardModelKWaveAdapter, VolumeCreationModelMode
    FieldOfViewCroppingProcessingComponent
from simpa.core.device_digital_twins import MSOTAcuityEcho

# FIXME temporary workaround for newest Intel architectures
import os
os.environ["KMP_DUPLICATE_LIB_OK"] = "TRUE"

VOLUME_TRANSDUCER_DIM_IN_MM = 75
VOLUME_PLANAR_DIM_IN_MM = 20

```

```

VOLUME_HEIGHT_IN_MM = 25
SPACING = 0.25
RANDOM_SEED = 4711

# TODO: Please make sure that a valid path_config.env file is located in your home directory
# point to the correct file in the PathManager().
path_manager = PathManager()

# If VISUALIZE is set to True, the simulation result will be plotted
VISUALIZE = True

def create_example_tissue():
    """
    This is a very simple example script of how to create a tissue definition.
    It contains a muscular background, an epidermis layer on top of the muscles
    and a blood vessel.
    """
    background_dictionary = Settings()
    background_dictionary[Tags.MOLECULE_COMPOSITION] = TISSUE_LIBRARY.constant(1e-10, 1e-10,
    background_dictionary[Tags.STRUCTURE_TYPE] = Tags.BACKGROUND

    muscle_dictionary = Settings()
    muscle_dictionary[Tags.PRIORITY] = 1
    muscle_dictionary[Tags.STRUCTURE_START_MM] = [0, 0, 0]
    muscle_dictionary[Tags.STRUCTURE_END_MM] = [0, 0, 100]
    muscle_dictionary[Tags.MOLECULE_COMPOSITION] = TISSUE_LIBRARY.muscle()
    muscle_dictionary[Tags.CONSIDER_PARTIAL_VOLUME] = True
    muscle_dictionary[Tags.ADHERE_TO_DEFORMATION] = True
    muscle_dictionary[Tags.STRUCTURE_TYPE] = Tags.HORIZONTAL_LAYER_STRUCTURE

    vessel_1_dictionary = Settings()
    vessel_1_dictionary[Tags.PRIORITY] = 3
    vessel_1_dictionary[Tags.STRUCTURE_START_MM] = [VOLUME_TRANSDUCER_DIM_IN_MM/2,
    0, 10]
    vessel_1_dictionary[Tags.STRUCTURE_END_MM] = [VOLUME_TRANSDUCER_DIM_IN_MM/2, VOLUME_PLAN
    vessel_1_dictionary[Tags.STRUCTURE_RADIUS_MM] = 3
    vessel_1_dictionary[Tags.MOLECULE_COMPOSITION] = TISSUE_LIBRARY.blood()
    vessel_1_dictionary[Tags.CONSIDER_PARTIAL_VOLUME] = True
    vessel_1_dictionary[Tags.STRUCTURE_TYPE] = Tags.CIRCULAR_TUBULAR_STRUCTURE

    epidermis_dictionary = Settings()
    epidermis_dictionary[Tags.PRIORITY] = 8
    epidermis_dictionary[Tags.STRUCTURE_START_MM] = [0, 0, 0]
    epidermis_dictionary[Tags.STRUCTURE_END_MM] = [0, 0, 0.5]
    epidermis_dictionary[Tags.MOLECULE_COMPOSITION] = TISSUE_LIBRARY.epidermis()
    epidermis_dictionary[Tags.CONSIDER_PARTIAL_VOLUME] = True
    epidermis_dictionary[Tags.ADHERE_TO_DEFORMATION] = True
    epidermis_dictionary[Tags.STRUCTURE_TYPE] = Tags.HORIZONTAL_LAYER_STRUCTURE

    tissue_dict = Settings()
    tissue_dict[Tags.BACKGROUND] = background_dictionary
    tissue_dict["muscle"] = muscle_dictionary
    tissue_dict["epidermis"] = epidermis_dictionary
    tissue_dict["vessel_1"] = vessel_1_dictionary
    return tissue_dict

# Seed the numpy random configuration prior to creating the global_settings file in
# order to ensure that the same volume

```

```

# is generated with the same random seed every time.

np.random.seed(RANDOM_SEED)
VOLUME_NAME = "CompletePipelineTestMSOT_"+str(RANDOM_SEED)

general_settings = {
    # These parameters set the general properties of the simulated volume
    Tags.RANDOM_SEED: RANDOM_SEED,
    Tags.VOLUME_NAME: "CompletePipelineTestMSOT_" + str(RANDOM_SEED),
    Tags.SIMULATION_PATH: path_manager.get_hdf5_file_save_path(),
    Tags.SPACING_MM: SPACING,
    Tags.DIM_VOLUME_Z_MM: VOLUME_HEIGHT_IN_MM,
    Tags.DIM_VOLUME_X_MM: VOLUME_TRANSDUCER_DIM_IN_MM,
    Tags.DIM_VOLUME_Y_MM: VOLUME_PLANAR_DIM_IN_MM,
    Tags.VOLUME_CREATOR: Tags.VOLUME_CREATOR_VERSATILE,
    Tags.GPU: True,

    # The following parameters set the optical forward model
    Tags.WAVELENGTHS: [700]
}
settings = Settings(general_settings)
np.random.seed(RANDOM_SEED)

settings.set_volume_creation_settings({
    Tags.STRUCTURES: create_example_tissue(),
    Tags.SIMULATE_DEFORMED_LAYERS: True
})

settings.set_optical_settings({
    Tags.OPTICAL_MODEL_NUMBER_PHOTONS: 1e7,
    Tags.OPTICAL_MODEL_BINARY_PATH: path_manager.get_mcx_binary_path(),
    Tags.ILLUMINATION_TYPE: Tags.ILLUMINATION_TYPE_MSOT_ACUITY_ECHO,
    Tags.LASER_PULSE_ENERGY_IN_MILLIJOULE: 50,
})

settings.set_acoustic_settings({
    Tags.ACOUSTIC_SIMULATION_3D: False,
    Tags.ACOUSTIC_MODEL_BINARY_PATH: path_manager.get_matlab_binary_path(),
    Tags.PROPERTY_ALPHA_POWER: 1.05,
    Tags.SENSOR_RECORD: "p",
    Tags.PMLInside: False,
    Tags.PMLSize: [31, 32],
    Tags.PMLAlpha: 1.5,
    Tags.PlotPML: False,
    Tags.RECORDMOVIE: False,
    Tags.MOVIE_NAME: "visualization_log",
    Tags.ACOUSTIC_LOG_SCALE: True
})

settings.set_reconstruction_settings({
    Tags.RECONSTRUCTION_PERFORM_BANDPASS_FILTERING: False,
    Tags.ACOUSTIC_MODEL_BINARY_PATH: path_manager.get_matlab_binary_path(),
    Tags.ACOUSTIC_SIMULATION_3D: True,
    Tags.PROPERTY_ALPHA_POWER: 1.05,
    Tags.TUKEY_WINDOW_ALPHA: 0.5,
    Tags.BANDPASS_CUTOFF_LOWPASS: int(8e6),
    Tags.BANDPASS_CUTOFF_HIGHPASS: int(0.1e6),
    Tags.RECONSTRUCTION_BMODE_AFTER_RECONSTRUCTION: True,
    Tags.RECONSTRUCTION_BMODE_METHOD: Tags.RECONSTRUCTION_BMODE_METHOD_HILBERT_TRANSFORM,
    Tags.RECONSTRUCTION_APODIZATION_METHOD: Tags.RECONSTRUCTION_APODIZATION_BOX,

```

```

    Tags.RECONSTRUCTION_MODE: Tags.RECONSTRUCTION_MODE_DIFFERENTIAL,
    Tags.SENSOR_RECORD: "p",
    Tags.PMLInside: False,
    Tags.PMLSize: [31, 32],
    Tags.PMLAlpha: 1.5,
    Tags.PlotPML: False,
    Tags.RECORDMOVIE: False,
    Tags.MOVIE_NAME: "visualization_log",
    Tags.ACOUSTIC_LOG_SCALE: True
})

settings["noise_initial_pressure"] = {
    Tags.NOISE_MEAN: 1,
    Tags.NOISE_STD: 0.01,
    Tags.NOISE_MODE: Tags.NOISE_MODE_MULTIPLICATIVE,
    Tags.DATA_FIELD: Tags.OPTICAL_MODEL_INITIAL_PRESSURE,
    Tags.NOISE_NON_NEGATIVITY_CONSTRAINT: True
}

settings["noise_time_series"] = {
    Tags.NOISE_STD: 1,
    Tags.NOISE_MODE: Tags.NOISE_MODE_ADDITIVE,
    Tags.DATA_FIELD: Tags.TIME_SERIES_DATA
}

device = MSOTAcuityEcho(device_position_mm=np.array([VOLUME_TRANSDUCER_DIM_IN_MM/2,
                                                    VOLUME_PLANAR_DIM_IN_MM/2,
                                                    0]))

device.update_settings_for_use_of_model_based_volume_creator(settings)

SIMULATION_PIPELINE = [
    VolumeCreationModelModelBasedAdapter(settings),
    OpticalForwardModelMcxAdapter(settings),
    GaussianNoiseProcessingComponent(settings, "noise_initial_pressure"),
    AcousticForwardModelKWaveAdapter(settings),
    GaussianNoiseProcessingComponent(settings, "noise_time_series"),
    FieldOfViewCroppingProcessingComponent(settings),
    ImageReconstructionModuleDelayAndSumAdapter(settings)
]

simulate(SIMULATION_PIPELINE, settings, device)

if Tags.WAVELENGTH in settings:
    WAVELENGTH = settings[Tags.WAVELENGTH]
else:
    WAVELENGTH = 700

if VISUALIZE:
    visualise_data(path_manager.get_hdf5_file_save_path() + "/" + VOLUME_NAME + ".hdf5", WAV
                  show_time_series_data=True,
                  show_initial_pressure=True,
                  show_absorption=False,
                  show_segmentation_map=False,
                  show_tissue_density=False,
                  show_reconstructed_data=True,
                  show_fluence=False,
                  log_scale=False)

```

## Defining custom tissue structures and properties

The file can be found in `simpa_examples/create_custom_tissues.py`:

```
"""
```

```
SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ
```

```
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI)
```

```
SPDX-License-Identifier: MIT
```

```
"""
```

```
from simpa.utils import MolecularCompositionGenerator
```

```
from simpa.utils import MOLECULE_LIBRARY
```

```
from simpa.utils import Molecule
```

```
from simpa.utils import Spectrum
```

```
from simpa.utils.libraries.spectra_library import ScatteringSpectrumLibrary, AnisotropySpectrumLibrary
```

```
import numpy as np
```

```
def create_custom_absorber():
```

```
    wavelengths = np.linspace(200, 1500, 100)
```

```
    absorber = Spectrum(spectrum_name="random absorber",
```

```
                        wavelengths=wavelengths,
```

```
                        values=np.random.random(
```

```
                            np.shape(wavelengths)))
```

```
    return absorber
```

```
def create_custom_chromophore(volume_fraction: float = 1.0):
```

```
    chromophore = Molecule(
```

```
        absorption_spectrum=create_custom_absorber(),
```

```
        volume_fraction=volume_fraction,
```

```
        scattering_spectrum=ScatteringSpectrumLibrary.CONSTANT_SCATTERING_ARBITRARY(40.0),
```

```
        anisotropy_spectrum=AnisotropySpectrumLibrary.CONSTANT_ANISOTROPY_ARBITRARY(0.9)
```

```
    )
```

```
    return chromophore
```

```
def create_custom_tissue_type():
```

```
    # First create an instance of a TissueSettingsGenerator
```

```
    tissue_settings_generator = MolecularCompositionGenerator()
```

```
    water_volume_fraction = 0.4
```

```
    blood_volume_fraction = 0.5
```

```
    custom_chromophore_volume_fraction = 0.1
```

```
    # The volume fraction within every tissue type should sum up to 1.
```

```
    oxygenation = 0.4
```

```
    # Then append chromophores that you want
```

```
    tissue_settings_generator.append(key="oxyhemoglobin",
```

```
                                    value=MOLECULE_LIBRARY.oxyhemoglobin(oxygenation * blood_volume_fraction))
```

```
    tissue_settings_generator.append(key="deoxyhemoglobin",
```

```
                                    value=MOLECULE_LIBRARY.deoxyhemoglobin((1 - oxygenation) * blood_volume_fraction))
```

```
    tissue_settings_generator.append(key="water",
```

```
                                    value=MOLECULE_LIBRARY.water(water_volume_fraction))
```

```
    tissue_settings_generator.append(key="custom",
```

```
                                    value=create_custom_chromophore(custom_chromophore_volume_fraction))
```

```
    return tissue_settings_generator
```



## Defining a custom digital device twin class

The file can be found in `simpa_examples/create_a_custom_digital_device_twin.py`:

```
"""
```

```
SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ
```

```
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI)
```

```
SPDX-License-Identifier: MIT
```

```
"""
```

```
# FIXME temporary workaround for newest Intel architectures
```

```
import os
```

```
os.environ["KMP_DUPLICATE_LIB_OK"] = "TRUE"
```

```
from simpa.core.device_digital_twins.digital_device_twin_base import PhotoacousticDevice
```

```
from simpa.core.device_digital_twins.devices.illumination_geometries.slit_illumination import
```

```
from simpa.core.device_digital_twins.devices.detection_geometries.linear_array import Linear
```

```
from simpa.utils import Settings, Tags
```

```
import numpy as np
```

```
class ExampleDeviceSlitIlluminationLinearDetector(PhotoacousticDevice):
```

```
    """
```

```
    This class represents a digital twin of a PA device with a slit as illumination next to
```

```
    """
```

```
    def __init__(self):
```

```
        super().__init__()
```

```
        self.set_detection_geometry(LinearArrayDetectionGeometry())
```

```
        self.add_illumination_geometry(SlitIlluminationGeometry())
```

```
if __name__ == "__main__":
```

```
    device = ExampleDeviceSlitIlluminationLinearDetector()
```

```
    settings = Settings()
```

```
    settings[Tags.DIM_VOLUME_X_MM] = 20
```

```
    settings[Tags.DIM_VOLUME_Y_MM] = 50
```

```
    settings[Tags.DIM_VOLUME_Z_MM] = 20
```

```
    settings[Tags.SPACING_MM] = 0.5
```

```
    settings[Tags.STRUCTURES] = {}
```

```
    x_dim = int(round(settings[Tags.DIM_VOLUME_X_MM]/settings[Tags.SPACING_MM]))
```

```
    z_dim = int(round(settings[Tags.DIM_VOLUME_Z_MM]/settings[Tags.SPACING_MM]))
```

```
    positions = device.get_detection_geometry().get_detector_element_positions_accounting_for
```

```
    detector_elements = device.get_detection_geometry().get_detector_element_orientations(gl
```

```
    positions = np.round(positions/settings[Tags.SPACING_MM]).astype(int)
```

```
    import matplotlib.pyplot as plt
```

```
    plt.scatter(positions[:, 0], positions[:, 2])
```

```
    plt.quiver(positions[:, 0], positions[:, 2], detector_elements[:, 0], detector_elements[
```

```
    plt.show()
```

## Defining custom tissue types

The file can be found in `simpa_examples/create_custom_tissues.py`:

```
"""
```

```
SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ
```

```
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI)
```

```
SPDX-License-Identifier: MIT
```

Load a segmentation mask and use it to simulate

```
"""
```

```
from simpa.utils import MolecularCompositionGenerator
from simpa.utils import MOLECULE_LIBRARY
from simpa.utils import Molecule
from simpa.utils import Spectrum
from simpa.utils.libraries.spectra_library import ScatteringSpectrumLibrary, AnisotropySpectrumLibrary
import numpy as np

def create_custom_absorber():
    wavelengths = np.linspace(200, 1500, 100)
    absorber = Spectrum(spectrum_name="random absorber",
                       wavelengths=wavelengths,
                       values=np.random.random(
                           np.shape(wavelengths)))

    return absorber

def create_custom_chromophore(volume_fraction: float = 1.0):
    chromophore = Molecule(
        absorption_spectrum=create_custom_absorber(),
        volume_fraction=volume_fraction,
        scattering_spectrum=ScatteringSpectrumLibrary.CONSTANT_SCATTERING_ARBITRARY(40.0),
        anisotropy_spectrum=AnisotropySpectrumLibrary.CONSTANT_ANISOTROPY_ARBITRARY(0.9)
    )
    return chromophore

def create_custom_tissue_type():
    # First create an instance of a TissueSettingsGenerator
    tissue_settings_generator = MolecularCompositionGenerator()

    water_volume_fraction = 0.4
    blood_volume_fraction = 0.5
    custom_chromophore_volume_fraction = 0.1
    # The volume fraction within every tissue type should sum up to 1.

    oxygenation = 0.4

    # Then append chromophores that you want
    tissue_settings_generator.append(key="oxyhemoglobin",
                                    value=MOLECULE_LIBRARY.oxyhemoglobin(oxygenation * blood_volume_fraction))
    tissue_settings_generator.append(key="deoxyhemoglobin",
                                    value=MOLECULE_LIBRARY.deoxyhemoglobin((1 - oxygenation) * blood_volume_fraction))
    tissue_settings_generator.append(key="water",
                                    value=MOLECULE_LIBRARY.water(water_volume_fraction))
    tissue_settings_generator.append(key="custom",
                                    value=create_custom_chromophore(custom_chromophore_volume_fraction))

    return tissue_settings_generator
```

## Load a segmentation mask and use it to simulate

The file can be found in `simpa_examples/segmentation_loader.py`:

```
"""
```

```
SPDX-FileCopyrightText: 2021 Computer Assisted Medical Interventions Group, DKFZ
SPDX-FileCopyrightText: 2021 VISION Lab, Cancer Research UK Cambridge Institute (CRUK CI)
```

*SPDX-License-Identifier: MIT*

"""

```

from simpa.core.simulation import simulate
from simpa.utils.settings import Settings
from simpa.utils import Tags, SegmentationClasses
import numpy as np
from skimage.data import shepp_logan_phantom
from simpa.utils.libraries.tissue_library import TISSUE_LIBRARY
from simpa.utils.libraries.molecule_library import MOLECULE_LIBRARY
from simpa.utils.libraries.tissue_library import MolecularCompositionGenerator
from simpa.visualisation.matplotlib_data_visualisation import visualise_data
from scipy.ndimage import zoom
from simpa.utils.path_manager import PathManager
from simpa.core.device_digital_twins import RSOMEExplorerP50
from simpa.core import VolumeCreationModuleSegmentationBasedAdapter, OpticalForwardModelMcxA

# FIXME temporary workaround for newest Intel architectures
import os
os.environ["KMP_DUPLICATE_LIB_OK"] = "TRUE"

# If VISUALIZE is set to True, the simulation result will be plotted
VISUALIZE = True

# TODO: Please make sure that a valid path_config.env file is located in your home directory
# point to the correct file in the PathManager().
path_manager = PathManager()

target_spacing = 1.0

label_mask = shepp_logan_phantom()

label_mask = np.digitize(label_mask, bins=np.linspace(0.0, 1.0, 11), right=True)

label_mask = np.reshape(label_mask, (400, 1, 400))

input_spacing = 0.2
segmentation_volume_tiled = np.tile(label_mask, (1, 128, 1))
segmentation_volume_mask = np.round(zoom(segmentation_volume_tiled, input_spacing/target_spacing,
                                          order=0)).astype(int)

def segmentation_class_mapping():
    ret_dict = dict()
    ret_dict[0] = TISSUE_LIBRARY.heavy_water()
    ret_dict[1] = TISSUE_LIBRARY.blood()
    ret_dict[2] = TISSUE_LIBRARY.epidermis()
    ret_dict[3] = TISSUE_LIBRARY.muscle()
    ret_dict[4] = TISSUE_LIBRARY.mediprene()
    ret_dict[5] = TISSUE_LIBRARY.ultrasound_gel()
    ret_dict[6] = TISSUE_LIBRARY.heavy_water()
    ret_dict[7] = (MolecularCompositionGenerator()
                  .append(MOLECULE_LIBRARY.oxyhemoglobin(0.01))
                  .append(MOLECULE_LIBRARY.deoxyhemoglobin(0.01))
                  .append(MOLECULE_LIBRARY.water(0.98))
                  .get_molecular_composition(SegmentationClasses.COUPLING_ARTIFACT))
    ret_dict[8] = TISSUE_LIBRARY.heavy_water()
    ret_dict[9] = TISSUE_LIBRARY.heavy_water()
    ret_dict[10] = TISSUE_LIBRARY.heavy_water()
    return ret_dict

```

```

settings = Settings()
settings[Tags.SIMULATION_PATH] = path_manager.get_hdf5_file_save_path()
settings[Tags.VOLUME_NAME] = "SegmentationTest"
settings[Tags.RANDOM_SEED] = 1234
settings[Tags.WAVELENGTHS] = [700]
settings[Tags.SPACING_MM] = target_spacing
settings[Tags.DIM_VOLUME_X_MM] = 400 / (target_spacing / input_spacing)
settings[Tags.DIM_VOLUME_Y_MM] = 128 / (target_spacing / input_spacing)
settings[Tags.DIM_VOLUME_Z_MM] = 400 / (target_spacing / input_spacing)

settings.set_volume_creation_settings({
    Tags.INPUT_SEGMENTATION_VOLUME: segmentation_volume_mask,
    Tags.SEGMENTATION_CLASS_MAPPING: segmentation_class_mapping(),
})

settings.set_optical_settings({
    Tags.OPTICAL_MODEL_NUMBER_PHOTONS: 1e7,
    Tags.OPTICAL_MODEL_BINARY_PATH: path_manager.get_mcx_binary_path(),
    Tags.ILLUMINATION_TYPE: Tags.ILLUMINATION_TYPE_MSOT_ACUITY_ECHO,
    Tags.LASER_PULSE_ENERGY_IN_MILLIJOULE: 50,
})

pipeline = [
    VolumeCreationModuleSegmentationBasedAdapter(settings),
    OpticalForwardModelMcxAdapter(settings)
]

simulate(pipeline, settings, RSOMEExplorerP50(element_spacing_mm=1.0))

if Tags.WAVELENGTH in settings:
    WAVELENGTH = settings[Tags.WAVELENGTH]
else:
    WAVELENGTH = 700

if VISUALIZE:
    visualise_data(path_manager.get_hdf5_file_save_path() + "/" + "SegmentationTest" + ".hdf5")

```

# Index

## A

ACOUSTIC\_LOG\_SCALE (simpa.core.devices.digital\_twins.devices.detection\_geometries.planar\_array.PlanarArrayDetectionGeometry attribute)  
(simpa.core.devices.digital\_twins.devices.illumination\_geometries.illumination\_geometry\_base.IlluminationGeometryBase attribute)  
ACOUSTIC\_MODEL (simpa.utils.tags.Tags attribute)  
ACOUSTIC\_MODEL\_BINARY\_PATH (simpa.utils.tags.Tags attribute)  
ACOUSTIC\_MODEL\_K\_WAVE (simpa.utils.tags.Tags attribute)  
ACOUSTIC\_MODEL\_OUTPUT\_NAME (simpa.utils.tags.Tags attribute)  
ACOUSTIC\_MODEL\_SETTINGS (simpa.utils.tags.Tags attribute)  
ACOUSTIC\_MODEL\_TEST (simpa.utils.tags.Tags attribute)  
ACOUSTIC\_SIMULATION\_3D (simpa.utils.tags.Tags attribute)  
AcousticForwardModelBaseAdapter (class in simpa.core.acoustic\_forward\_module)  
AcousticForwardModelKWaveAdapter (class in simpa.core.acoustic\_forward\_module.acoustic\_forward\_model\_k\_wave\_adapter)  
ADHERE\_TO\_DEFORMATION (simpa.core.volume\_creation\_module.VolumeCreatorModuleBase attribute)  
apply\_b\_mode() (in module simpa.core.reconstruction\_module.reconstruction\_utils)

## B

Background (class in simpa.utils.libraries.structure\_library)  
BACKGROUND (simpa.utils.tags.Tags attribute)  
BANDPASS\_CUTOFF\_HIGHPASS (simpa.utils.tags.Tags attribute)  
BANDPASS\_CUTOFF\_LOWPASS (simpa.utils.tags.Tags attribute)  
bandpass\_filtering() (in module simpa.core.reconstruction\_module.reconstruction\_utils)  
blood() (simpa.utils.libraries.tissue\_library.TissueLibrary method)  
bone() (simpa.utils.libraries.tissue\_library.TissueLibrary method)

## C

calculate\_gruneisen\_parameter\_from\_temperature() (in module simpa.utils.calculate)  
calculate\_oxygenation() (in module simpa.utils.calculate)

check\_settings\_prerequisites() (simpa.core.devices.digital\_twins.devices.detection\_geometries.curved\_array.CurvedArrayDetectionGeometry method)

(simpa.core.devices.digital\_twins.devices.detection\_geometries.linear\_array.LinearArrayDetectionGeometry method)  
CIRCULAR\_TUBULAR\_STRUCTURE (simpa.utils.tags.Tags attribute)  
CircularTubularStructure (class in simpa.utils.libraries.structure\_library)  
CONSIDER\_PARTIAL\_VOLUME (simpa.utils.tags.Tags attribute)  
constant() (simpa.utils.libraries.tissue\_library.TissueLibrary method)  
create\_deformation\_settings() (in module simpa.utils.deformation\_manager)  
create\_simulation\_volume() (simpa.core.volume\_creation\_module.volume\_creation\_module\_model\_based\_adapter.VolumeCreationModelModelBasedAdapter method)

(simpa.core.volume\_creation\_module.volume\_creation\_module\_segmentation\_based\_adapter.VolumeCreationModelModelSegmentationBasedAdapter method)  
create\_spline\_for\_range() (in module simpa.utils.calculate)  
critical() (simpa.log.file\_logger.Logger method)  
CurvedArrayDetectionGeometry (class in simpa.core.devices.digital\_twins.devices.detection\_geometries.curved\_array)

## D

DATA\_FIELD (simpa.utils.tags.Tags attribute)  
debug() (simpa.log.file\_logger.Logger method)  
default() (simpa.io\_handling.serialization.SIMPAJSONSerializer method)  
DEFORMATION\_X\_COORDINATES\_MM (simpa.utils.tags.Tags attribute)  
DEFORMATION\_Y\_COORDINATES\_MM (simpa.utils.tags.Tags attribute)  
DEFORMATION\_Z\_ELEVATIONS\_MM (simpa.utils.tags.Tags attribute)  
DEFORMED\_LAYERS\_SETTINGS (simpa.utils.tags.Tags attribute)  
dermis() (simpa.utils.libraries.tissue\_library.TissueLibrary method)

`detect_local_path_config()`  
(`simpa.utils.path_manager.PathManager` method)

`DetectionGeometryBase` (class in `simpa.core.device_digital_twins.devices.detection_geometries.detection_geometry_base`)

`DETECTOR_ELEMENT_WIDTH_MM`  
(`simpa.utils.tags.Tags` attribute)

`DIGITAL_DEVICE` (`simpa.utils.tags.Tags` attribute)

`DIGITAL_DEVICE_MSOT_ACUITY`  
(`simpa.utils.tags.Tags` attribute)

`DIGITAL_DEVICE_MSOT_INVISION`  
(`simpa.utils.tags.Tags` attribute)

`DIGITAL_DEVICE_POSITION` (`simpa.utils.tags.Tags` attribute)

`DIGITAL_DEVICE_RSOM` (`simpa.utils.tags.Tags` attribute)

`DIGITAL_DEVICE_SLIT_ILLUMINATION_LINEAR_DETECTOR` (`simpa.utils.tags.Tags` attribute)

`DIM_VOLUME_X_MM` (`simpa.utils.tags.Tags` attribute)

`DIM_VOLUME_Y_MM` (`simpa.utils.tags.Tags` attribute)

`DIM_VOLUME_Z_MM` (`simpa.utils.tags.Tags` attribute)

`DOWNSCALE_FACTOR` (`simpa.utils.tags.Tags` attribute)

## E

`ELLIPTICAL_TUBULAR_STRUCTURE` (`simpa.core.device_digital_twins.devices.detection_geometries.linear_array.LinearArrayDetectionGeometryBase` method)

`EllipticalTubularStructure` (class in `simpa.utils.libraries.structure_library`)

`epidermis()`  
(`simpa.utils.libraries.tissue_library.TissueLibrary` method)

`EPS` (in module `simpa.utils.constants`)

`error()` (`simpa.log.file_logger.Logger` method)

## F

`fill_internal_volume()` (`simpa.utils.libraries.structure_library.GeometricalStructure` method)

`fill_vessel_volume()` (`simpa.utils.libraries.structure_library.VesselStructure` method)

`forward_model()` (`simpa.core.acoustic_forward_module.acoustic_forward_module_k_wave_adapter.AcousticForwardModelKWWaveAdapter` method)

`forward_model()` (`simpa.core.acoustic_forward_module.acoustic_forward_module_base_adapter.AcousticForwardModelBaseAdapter` method)

`forward_model()` (`simpa.core.acoustic_forward_module.acoustic_forward_module_mcx_adapter.OpticalForwardModelMCXAdapter` method)

`forward_model()` (`simpa.core.acoustic_forward_module.acoustic_forward_module_base_adapter.AcousticForwardModelBaseAdapter` method)

## G

`GammaNoiseProcessingComponent` (class in `simpa.core.processing_components.noise.gamma_noise`)

`GaussianNoiseProcessingComponent` (class in `simpa.core.processing_components.noise.gaussian_noise`)

`generate_dict_path()` (in module `simpa.utils.dict_path_manager`)

`GeometricalStructure` (class in `simpa.utils.libraries.structure_library`)

`get_acoustic_properties()` (`simpa.core.reconstruction_module.reconstruction_module_time_reversal_adapter.ReconstructionModuleTimeReversalAdapter` method)

`get_acoustic_settings()` (`simpa.utils.settings.Settings` method)

`get_apodization_factor()` (in module `simpa.core.reconstruction_module.reconstruction_utils`)

`get_blood_volume_fractions()`  
(`simpa.utils.libraries.tissue_library.TissueLibrary` method)

`get_detector_element_orientations()` (`simpa.core.device_digital_twins.devices.detection_geometries.curved_array.CurvedArrayDetectionGeometry` method)

`get_detector_element_orientations()` (`simpa.core.device_digital_twins.devices.detection_geometries.detection_geometry_base.DetectionGeometryBase` method)

`get_detector_element_orientations()` (`simpa.core.device_digital_twins.devices.detection_geometries.linear_array.LinearArrayDetectionGeometryBase` method)

`get_detector_element_orientations()` (`simpa.core.device_digital_twins.devices.detection_geometries.planar_array.PlanarArrayDetectionGeometryBase` method)

`get_detector_element_positions_accounting_for_device_position_mm()` (`simpa.core.device_digital_twins.devices.detection_geometries.detection_geometry_base.DetectionGeometryBase` method)

`get_detector_element_positions_accounting_for_field_of_view()` (`simpa.core.device_digital_twins.devices.detection_geometries.detection_geometry_base.DetectionGeometryBase` method)

`get_detector_element_positions_base_mm()` (`simpa.core.device_digital_twins.devices.detection_geometries.curved_array.CurvedArrayDetectionGeometry` method)

`get_detector_element_positions_base_mm()` (`simpa.core.device_digital_twins.devices.detection_geometries.detection_geometry_base.DetectionGeometryBase` method)

`get_detector_element_positions_base_mm()` (`simpa.core.device_digital_twins.devices.detection_geometries.linear_array.LinearArrayDetectionGeometryBase` method)

`get_detector_element_positions_base_mm()` (`simpa.core.device_digital_twins.devices.detection_geometries.planar_array.PlanarArrayDetectionGeometryBase` method)

`get_enclosed_indices()`

(`simpa.utils.libraries.structure_library.BackgroundForwardModelMCXAdapter` method)

(`simpa.utils.libraries.structure_library.CircularTubularStructure` method)

(`simpa.utils.libraries.structure_library.EllipticalTubularStructure` method)



(simpa.utils.libraries.structure_library.GeometricalStructure method)	(simpa.utils.libraries.structure_library.VesselStructure method)
(simpa.utils.libraries.structure_library.HorizontalLayerStructure method)	get_reconstruction_settings() (simpa.utils.settings.Settings method)
(simpa.utils.libraries.structure_library.ParallelepipedStructure method)	get_value_for_wavelength() (simpa.utils.libraries.spectra_library.Spectrum method)
(simpa.utils.libraries.structure_library.RectangularCuboidStructure method)	get_value_over_wavelength() (simpa.utils.libraries.spectra_library.Spectrum method)
(simpa.utils.libraries.structure_library.SphericalStructure method)	get_volume_creation_settings() (simpa.utils.settings.Settings method)
(simpa.utils.libraries.structure_library.VesselStructure method)	GPU (simpa.utils.tags.Tags attribute)
get_field_of_view_extent_mm() (simpa.core.device_digital_twins.devices.detection_geometries.curved_array.CurvedArrayDetectionGeometry method)	
ital_twins.devices.detection_geometries.linear_array.LinearArrayDetectionGeometry method)	<b>H</b>
ital_twins.devices.detection_geometries.planar_array.PlanarArrayDetectionGeometry method)	<b>HORIZONTAL_LAYER_STRUCTURE</b> (simpa.utils.tags.Tags attribute)
ital_twins.devices.illumination_geometries.pencil_array_illumination.PencilArrayIlluminationGeometry method)	HorizontalLayerStructure (class in simpa.utils.libraries.structure_library)
ital_twins.devices.illumination_geometries.pencil_beam_illumination.PencilBeamIlluminationGeometry method)	<b>I</b>
ital_twins.devices.illumination_geometries.slit_illumination.SlitIlluminationGeometry method)	ILLUMINATION_DIRECTION (simpa.utils.tags.Tags attribute)
ital_twins.devices.pa_devices.ithera_msot_invision.InVision256TF method)	ILLUMINATION_PARAM1 (simpa.utils.tags.Tags attribute)
get_functional_from_deformation_settings() (in module simpa.utils.deformation_manager)	ILLUMINATION_PARAM2 (simpa.utils.tags.Tags attribute)
get_mcx_illuminator_definition() (simpa.core.device_digital_twins.devices.illumination_geometries.illumination_geometry_base.IlluminationGeometryBase method)	ILLUMINATION_POSITION (simpa.utils.tags.Tags attribute)
ital_twins.devices.illumination_geometries.pencil_array_illumination.PencilArrayIlluminationGeometry method)	ILLUMINATION_TYPE (simpa.utils.tags.Tags attribute)
ital_twins.devices.illumination_geometries.pencil_beam_illumination.PencilBeamIlluminationGeometry method)	ILLUMINATION_TYPE_DISK (simpa.utils.tags.Tags attribute)
ital_twins.devices.illumination_geometries.slit_illumination.SlitIlluminationGeometry method)	ILLUMINATION_TYPE_DKFZ_PAUS (simpa.utils.tags.Tags attribute)
get_optical_settings() (simpa.utils.settings.Settings method)	ILLUMINATION_TYPE_FOURIER (simpa.utils.tags.Tags attribute)
get_params_from_settings() (simpa.utils.libraries.structure_library.Background method)	ILLUMINATION_TYPE_FOURIER_X (simpa.utils.tags.Tags attribute)
(simpa.utils.libraries.structure_library.CircularTubularStructure method)	ILLUMINATION_TYPE_FOURIER_X_2D (simpa.utils.tags.Tags attribute)
(simpa.utils.libraries.structure_library.EllipticalTubularStructure method)	ILLUMINATION_TYPE_GAUSSIAN (simpa.utils.tags.Tags attribute)
(simpa.utils.libraries.structure_library.GeometricalStructure method)	ILLUMINATION_TYPE_IPASC_DEFINITION (simpa.utils.tags.Tags attribute)
(simpa.utils.libraries.structure_library.HorizontalLayerStructure method)	ILLUMINATION_TYPE_MSOT_ACUITY_ECHO (simpa.utils.tags.Tags attribute)
(simpa.utils.libraries.structure_library.ParallelepipedStructure method)	ILLUMINATION_TYPE_MSOT_INVISION (simpa.utils.tags.Tags attribute)
(simpa.utils.libraries.structure_library.RectangularCuboidStructure method)	ILLUMINATION_TYPE_PATTERN (simpa.utils.tags.Tags attribute)
(simpa.utils.libraries.structure_library.SphericalStructure method)	ILLUMINATION_TYPE_PATTERN_3D (simpa.utils.tags.Tags attribute)

ILLUMINATION\_TYPE\_PENCIL (simpa.utils.tags.Tags attribute)

ILLUMINATION\_TYPE\_PENCILARRAY (simpa.utils.tags.Tags attribute)

ILLUMINATION\_TYPE\_PLANAR (simpa.utils.tags.Tags attribute)

ILLUMINATION\_TYPE\_RING (simpa.utils.tags.Tags attribute)

ILLUMINATION\_TYPE\_SLIT (simpa.utils.tags.Tags attribute)

IlluminationGeometryBase (class in simpa.core.device\_digital\_twins.devices.illumination\_geometries.illumination\_geometry\_base)

IMAGE\_PROCESSING (simpa.utils.tags.Tags attribute)

ImageReconstructionModuleDelayAndSumAdapter (class in simpa.core.reconstruction\_module.reconstruction\_module\_delay\_and\_sum\_adapter)

ImageReconstructionModuleDelayMultiplyAndSumAdapter (class in simpa.core.reconstruction\_module.reconstruction\_module\_delay\_multiply\_and\_sum\_adapter)

ImageReconstructionModuleSignedDelayMultiplyAndSumAdapter (class in simpa.core.reconstruction\_module.reconstruction\_module\_signed\_delay\_multiply\_and\_sum\_adapter)

info() (simpa.log.file\_logger.Logger method)

INPUT\_SEGMENTATION\_VOLUME (simpa.utils.tags.Tags attribute)

InVision256TF (class in simpa.core.device\_digital\_twins.devices.pa\_devices.ithera\_msot\_invision)

ITERATIVE\_qPAI\_RESULT (simpa.utils.tags.Tags attribute)

ITERATIVE\_RECONSTRUCTION\_CONSTANT\_REGULARIZATION (simpa.utils.tags.Tags attribute)

ITERATIVE\_RECONSTRUCTION\_MAX\_ITERATION\_NUMBER (simpa.utils.tags.Tags attribute)

ITERATIVE\_RECONSTRUCTION\_REGULARIZATION\_SIGMA (simpa.utils.tags.Tags attribute)

ITERATIVE\_RECONSTRUCTION\_SAVE\_INTERMEDIATE\_RESULTS (simpa.utils.tags.Tags attribute)

ITERATIVE\_RECONSTRUCTION\_STOPPING\_LEVEL (simpa.utils.tags.Tags attribute)

## K

K\_WAVE\_SPECIFIC\_DT (simpa.utils.tags.Tags attribute)

K\_WAVE\_SPECIFIC\_NT (simpa.utils.tags.Tags attribute)

## L

LASER\_PULSE\_ENERGY\_IN\_MILLIJOULE (simpa.utils.tags.Tags attribute)

LINEAR\_UNMIXING\_RESULT (simpa.utils.tags.Tags attribute)

LinearArrayDetectionGeometry (class in simpa.core.device\_digital\_twins.devices.detection\_geometries.linear\_array)

LOAD\_AND\_SAVE\_HDF5\_FILE\_AT\_THE\_END\_OF\_SIMULATION\_TO\_MINIMISE\_FILESIZE (simpa.utils.tags.Tags attribute)

load\_hdf5() (in module simpa.io\_handling.io\_hdf5)

Logger (class in simpa.log.file\_logger)

## M

MAX\_DEFORMATION\_MM (simpa.utils.tags.Tags attribute)

MCX\_ASSUMED\_ANISOTROPY (simpa.utils.tags.Tags attribute)

MCX\_SEED (simpa.utils.tags.Tags attribute)

MEDIUM\_TEMPERATURE\_CELCIUS (simpa.utils.tags.Tags attribute)

min\_max\_normalization() (in module simpa.utils.calculate)

MolecularComposition (class in simpa.utils.libraries.molecule\_library)

MolecularCompositionGenerator (class in simpa.utils.libraries.tissue\_library)

MOLECULE\_COMPOSITION (simpa.utils.tags.Tags attribute)

MorphologicalTissueProperties (class in simpa.utils.libraries.literature\_values)

MOVIE\_NAME (simpa.utils.tags.Tags attribute)

MSOTAcuityEcho (class in simpa.core.device\_digital\_twins.devices.pa\_devices.ithera\_msot\_acuity)

muscle() (simpa.utils.libraries.tissue\_library.TissueLibrary method)

## N

NOISE\_FREQUENCY (simpa.utils.tags.Tags attribute)

NOISE\_MAX (simpa.utils.tags.Tags attribute)

NOISE\_MEAN (simpa.utils.tags.Tags attribute)

NOISE\_MIN (simpa.utils.tags.Tags attribute)

NOISE\_MODE (simpa.utils.tags.Tags attribute)

NOISE\_MODE\_ADDITIVE (simpa.utils.tags.Tags attribute)

NOISE\_MODE\_MULTIPLICATIVE (simpa.utils.tags.Tags attribute)



NOISE\_NON\_NEGATIVITY\_CONSTRAINT  
(simpa.utils.tags.Tags attribute)

NOISE\_SCALE (simpa.utils.tags.Tags attribute)

NOISE\_SHAPE (simpa.utils.tags.Tags attribute)

NOISE\_STD (simpa.utils.tags.Tags attribute)

## O

OPTICAL\_MODEL (simpa.utils.tags.Tags attribute)

OPTICAL\_MODEL\_BINARY\_PATH  
(simpa.utils.tags.Tags attribute)

OPTICAL\_MODEL\_FLUENCE (simpa.utils.tags.Tags attribute)

OPTICAL\_MODEL\_ILLUMINATION\_GEOMETRY\_JSON\_FILE (simpa.utils.tags.Tags attribute)

OPTICAL\_MODEL\_INITIAL\_PRESSURE  
(simpa.utils.tags.Tags attribute)

OPTICAL\_MODEL\_MCX (simpa.utils.tags.Tags attribute)

OPTICAL\_MODEL\_NUMBER\_PHOTONS  
(simpa.utils.tags.Tags attribute)

OPTICAL\_MODEL\_OUTPUT\_NAME  
(simpa.utils.tags.Tags attribute)

OPTICAL\_MODEL\_SETTINGS (simpa.utils.tags.Tags attribute)

OPTICAL\_MODEL\_TEST (simpa.utils.tags.Tags attribute)

OPTICAL\_MODEL\_UNITS (simpa.utils.tags.Tags attribute)

OpticalForwardModelMcxAdapter (class in simpa.core.optical\_simulation\_module.optical\_forward\_model\_mcx\_adapter)

OpticalForwardModuleBase (class in simpa.core.optical\_simulation\_module)

OpticalTissueProperties (class in simpa.utils.libraries.literature\_values)

ORIGINAL\_DATA (simpa.utils.tags.Tags attribute)

## P

PARALLELEPIPED\_STRUCTURE  
(simpa.utils.tags.Tags attribute)

ParallelepipedStructure (class in simpa.utils.libraries.structure\_library)

PathManager (class in simpa.utils.path\_manager)

PencilArrayIlluminationGeometry (class in simpa.core.device\_digital\_twins.devices.illumination\_geometries.pencil\_array\_illumination)

PencilBeamIlluminationGeometry (class in simpa.core.device\_digital\_twins.devices.illumination\_geometries.pencil\_beam\_illumination)

PlanarArrayDetectionGeometry (class in simpa.core.device\_digital\_twins.devices.detection\_geometries.planar\_array)

PlotPML (simpa.utils.tags.Tags attribute)

PMLAlpha (simpa.utils.tags.Tags attribute)

PMLInside (simpa.utils.tags.Tags attribute)

PMLSize (simpa.utils.tags.Tags attribute)

PoissonNoiseProcessingComponent (class in simpa.core.processing\_components.noise.poisson\_noise)

PRIORITY (simpa.utils.tags.Tags attribute)

ProcessingComponent (class in simpa.core.processing\_components)

properties\_for\_wavelength() (simpa.utils.libraries.structure\_library.GeometricalStructure method)

PROPERTY\_ABSORPTION\_PER\_CM  
(simpa.utils.tags.Tags attribute)

PROPERTY\_ALPHA\_COEFF (simpa.utils.tags.Tags attribute)

PROPERTY\_ALPHA\_POWER (simpa.utils.tags.Tags attribute)

PROPERTY\_ANISOTROPY (simpa.utils.tags.Tags attribute)

PROPERTY\_DENSITY (simpa.utils.tags.Tags attribute)

PROPERTY\_DIRECTIVITY\_ANGLE  
(simpa.utils.tags.Tags attribute)

PROPERTY\_GRUNEISEN\_PARAMETER  
(simpa.utils.tags.Tags attribute)

PROPERTY\_INTRINSIC\_EULER\_ANGLE  
(simpa.utils.tags.Tags attribute)

PROPERTY\_OXYGENATION (simpa.utils.tags.Tags attribute)

PROPERTY\_SCATTERING\_PER\_CM  
(simpa.utils.tags.Tags attribute)

PROPERTY\_SEGMENTATION (simpa.utils.tags.Tags attribute)

PROPERTY\_SENSOR\_MASK (simpa.utils.tags.Tags attribute)

PROPERTY\_SPEED\_OF\_SOUND  
(simpa.utils.tags.Tags attribute)

## R

RANDOM\_SEED (simpa.utils.tags.Tags attribute)

randomize\_uniform() (in module simpa.utils.calculate)

reconstruct\_delay\_and\_sum\_pytorch() (in module simpa.core.reconstruction\_module.reconstruction\_module\_delay\_and\_sum\_adapter)

reconstruct\_delay\_multiply\_and\_sum\_pytorch() (in module simpa.core.reconstruction\_module.reconstruction\_module\_delay\_multiply\_and\_sum\_adapter)

reconstruct_signed_delay_multiply_and_sum_pytorch() (in module simpa.core.reconstruction_module.reconstruction_module_signed_delay_multiply_and_sum_adapter)	RECONSTRUCTION_MITK_BINARY_PATH (simpa.utils.tags.Tags attribute)
RECONSTRUCTED_DATA (simpa.utils.tags.Tags attribute)	RECONSTRUCTION_MITK_SETTINGS_XML (simpa.utils.tags.Tags attribute)
RECONSTRUCTED_DATA_NOISE (simpa.utils.tags.Tags attribute)	RECONSTRUCTION_MODE (simpa.utils.tags.Tags attribute)
RECONSTRUCTION_ALGORITHM (simpa.utils.tags.Tags attribute)	RECONSTRUCTION_MODE_DIFFERENTIAL (simpa.utils.tags.Tags attribute)
reconstruction_algorithm() (simpa.core.reconstruction_module.reconstruction_module_delay_and_sum_adapter.ImageReconstructionModuleDelayAndSumAdapter method)	RECONSTRUCTION_MODE_FULL (simpa.utils.tags.Tags attribute)
reconstruction_module_delay_multiply_and_sum_adapter.ImageReconstructionModuleDelayMultiplyAndSumAdapter method)	RECONSTRUCTION_MODE_PRESSURE (simpa.utils.tags.Tags attribute)
reconstruction_module_signed_delay_multiply_and_sum_adapter.ImageReconstructionModuleSignedDelayMultiplyAndSumAdapter (simpa.utils.tags.Tags attribute)	reconstruction_mode_transformation() (in module simpa.core.reconstruction_module.reconstruction_utils)
reconstruction_module_time_reversal_adapter.ReconstructionModuleTimeReversalAdapter method)	RECONSTRUCTION_MODEL_SETTINGS (simpa.utils.tags.Tags attribute)
ReconstructionAdapterBase method)	RECONSTRUCTION_OUTPUT_NAME (simpa.utils.tags.Tags attribute)
RECONSTRUCTION_ALGORITHM_DAS (simpa.utils.tags.Tags attribute)	RECONSTRUCTION_PERFORM_BANDPASS_FILTERING (simpa.utils.tags.Tags attribute)
RECONSTRUCTION_ALGORITHM_DMAS (simpa.utils.tags.Tags attribute)	ReconstructionAdapterBase (class in simpa.core.reconstruction_module)
RECONSTRUCTION_ALGORITHM_PYTORCH_DAS (simpa.utils.tags.Tags attribute)	ReconstructionModuleTimeReversalAdapter (class in simpa.core.reconstruction_module.reconstruction_module_time_reversal_adapter)
RECONSTRUCTION_ALGORITHM_SDMAS (simpa.utils.tags.Tags attribute)	RECORDMOVIE (simpa.utils.tags.Tags attribute)
RECONSTRUCTION_ALGORITHM_TEST (simpa.utils.tags.Tags attribute)	RECTANGULAR_CUBOID_STRUCTURE (simpa.utils.tags.Tags attribute)
RECONSTRUCTION_ALGORITHM_TIME_REVERSAL (simpa.utils.tags.Tags attribute)	RectangularCuboidStructure (class in simpa.utils.libraries.structure_library)
RECONSTRUCTION_APODIZATION_BOX (simpa.utils.tags.Tags attribute)	reorder_time_series_data() (simpa.core.reconstruction_module.reconstruction_module_time_reversal_adapter.ReconstructionModuleTimeReversalAdapter method)
RECONSTRUCTION_APODIZATION_HAMMING (simpa.utils.tags.Tags attribute)	rotation() (in module simpa.utils.calculate)
RECONSTRUCTION_APODIZATION_HANN (simpa.utils.tags.Tags attribute)	rotation_matrix_between_vectors() (in module simpa.utils.calculate)
RECONSTRUCTION_APODIZATION_METHOD (simpa.utils.tags.Tags attribute)	rotation_x() (in module simpa.utils.calculate)
RECONSTRUCTION_BMODE_AFTER_RECONSTRUCTION (simpa.utils.tags.Tags attribute)	rotation_y() (in module simpa.utils.calculate)
RECONSTRUCTION_BMODE_BEFORE_RECONSTRUCTION (simpa.utils.tags.Tags attribute)	rotation_z() (in module simpa.utils.calculate)
RECONSTRUCTION_BMODE_METHOD (simpa.utils.tags.Tags attribute)	RSOMExplorerP50 (class in simpa.core.device_digital_twins.devices.pa_devices.ithera_rsom)
RECONSTRUCTION_BMODE_METHOD_ABS (simpa.core.processing_components.noise.gamma_noise.GammaNoiseProcessing method)	run() (simpa.core.acoustic_forward_module.AcousticForwardModelBaseAdapter method)
RECONSTRUCTION_BMODE_METHOD_HILBERT_TRANSFORM (simpa.core.processing_components.noise.gaussian_noise.GaussianNoiseProcessing method)	
RECONSTRUCTION_INVERSE_CRIME (simpa.utils.tags.Tags attribute)	

processing_components.noise.poisson_noise.PoissonNoiseProcessingComponentBase method)	set_optical_settings() (simpa.utils.settings.Settings method)
processing_components.noise.salt_and_pepper_noise.SaltAndPepperNoiseProcessingComponentBase method)	Settings (class in simpa.utils.settings)
processing_components.noise.uniform_noise.UniformNoiseProcessingComponentBase method)	SETTINGS (class in simpa.utils.tags.Tags attribute)
reconstruction_module.ReconstructionAdapterBase method)	SETTINGS_JSON (simpa.utils.tags.Tags attribute)
volume_creation_module.VolumeCreatorModuleBase method)	SETTINGS_JSON_PATH (simpa.utils.tags.Tags attribute)
<b>S</b>	
SaltAndPepperNoiseProcessingComponent (class in simpa.core.processing_components.noise.salt_and_pepper_noise)	simpa.core.acoustic_forward_module (module)
save_hdf5() (in module simpa.io_handling.io_hdf5)	simpa.core.acoustic_forward_module.acoustic_forward_module_k_wave_adapter (module)
SaveFilePaths (class in simpa.utils.constants)	simpa.core.device_digital_twins.devices.detection_geometries.curved_array (module)
SEGMENTATION_CLASS_MAPPING (simpa.utils.tags.Tags attribute)	simpa.core.device_digital_twins.devices.detection_geometries.detection_geometry_base (module)
SegmentationClasses (class in simpa.utils.constants)	simpa.core.device_digital_twins.devices.detection_geometries.linear_array (module)
SENSOR_BANDWIDTH_PERCENT (simpa.utils.tags.Tags attribute)	simpa.core.device_digital_twins.devices.detection_geometries.planar_array (module)
SENSOR_CENTER_FREQUENCY_HZ (simpa.utils.tags.Tags attribute)	simpa.core.device_digital_twins.devices.illumination_geometries.illumination_geometry_base (module)
SENSOR_CONCAVE (simpa.utils.tags.Tags attribute)	simpa.core.device_digital_twins.devices.illumination_geometries.pencil_array_illumination (module)
SENSOR_DIRECTIVITY_PATTERN (simpa.utils.tags.Tags attribute)	simpa.core.device_digital_twins.devices.illumination_geometries.pencil_beam_illumination (module)
SENSOR_DIRECTIVITY_SIZE_M (simpa.utils.tags.Tags attribute)	simpa.core.device_digital_twins.devices.illumination_geometries.slit_illumination (module)
SENSOR_ELEMENT_POSITIONS (simpa.utils.tags.Tags attribute)	simpa.core.device_digital_twins.devices.pa_devices.ithera_msot_acuity (module)
SENSOR_LINEAR (simpa.utils.tags.Tags attribute)	simpa.core.device_digital_twins.devices.pa_devices.ithera_msot_invision (module)
SENSOR_NUM_ELEMENTS (simpa.utils.tags.Tags attribute)	simpa.core.device_digital_twins.devices.pa_devices.ithera_rsom (module)
SENSOR_NUM_USED_ELEMENTS (simpa.utils.tags.Tags attribute)	simpa.core.optical_simulation_module (module)
SENSOR_PITCH_MM (simpa.utils.tags.Tags attribute)	simpa.core.optical_simulation_module.optical_forward_model_mcx_adapter (module)
SENSOR_RADIUS_MM (simpa.utils.tags.Tags attribute)	simpa.core.processing_components (module)
SENSOR_RECORD (simpa.utils.tags.Tags attribute)	simpa.core.processing_components.noise.gamma_noise (module)
SENSOR_SAMPLING_RATE_MHZ (simpa.utils.tags.Tags attribute)	simpa.core.processing_components.noise.gaussian_noise (module)
serialize() (simpa.io_handling.serialization.SIMPASerializer method)	simpa.core.processing_components.noise.poisson_noise (module)
set_acoustic_settings() (simpa.utils.settings.Settings method)	simpa.core.processing_components.noise.salt_and_pepper_noise (module)
set_optical_settings() (simpa.utils.settings.Settings method)	simpa.core.processing_components.noise.uniform_noise (module)
set_reconstruction_settings() (simpa.utils.settings.Settings method)	simpa.core.reconstruction_module (module)
	simpa.core.reconstruction_module.reconstruction_module_delay_and_sum_adapter (module)

<code>simpa.core.reconstruction_module.reconstruction_module_delay_multiply_and_sum_adapter</code> (module)	<code>SIMULATION_PROPERTIES</code> (simpa.utils.tags.Tags attribute)
<code>simpa.core.reconstruction_module.reconstruction_module_signed_delay_multiply_and_sum_adapter</code> (module)	<code>SIMULATIONS</code> (simpa.utils.tags.Tags attribute)
<code>simpa.core.reconstruction_module.reconstruction_module_time_reversal_adapter</code> (module)	<code>SlitIlluminationGeometry</code> (class in <code>simpa.core.device_digital_twins.devices.illumination_geometries.slit_illumination</code> )
<code>simpa.core.reconstruction_module.reconstruction_utils</code> (module)	<code>SPACING_MM</code> (simpa.utils.tags.Tags attribute)
<code>simpa.core.simulation</code> (module)	<code>Spectrum</code> (class in <code>simpa.utils.libraries.spectra_library</code> )
<code>simpa.core.volume_creation_module</code> (module)	<code>SPHERICAL_STRUCTURE</code> (simpa.utils.tags.Tags attribute)
<code>simpa.core.volume_creation_module.volume_creation_module_model_based_adapter</code> (module)	<code>SphericalStructure</code> (class in <code>simpa.utils.libraries.structure_library</code> )
<code>simpa.core.volume_creation_module.volume_creation_module_segmentation_based_adapter</code> (module)	<code>StandardProperties</code> (class in <code>simpa.utils.libraries.literature_values</code> )
<code>simpa.io_handling</code> (module)	<code>STRUCTURE_BIFURCATION_LENGTH_MM</code> (simpa.utils.tags.Tags attribute)
<code>simpa.io_handling.io_hdf5</code> (module)	<code>STRUCTURE_CURVATURE_FACTOR</code> (simpa.utils.tags.Tags attribute)
<code>simpa.io_handling.serialization</code> (module)	<code>STRUCTURE_DIRECTION</code> (simpa.utils.tags.Tags attribute)
<code>simpa.log</code> (module)	<code>STRUCTURE_ECCENTRICITY</code> (simpa.utils.tags.Tags attribute)
<code>simpa.log.file_logger</code> (module)	<code>STRUCTURE_END_MM</code> (simpa.utils.tags.Tags attribute)
<code>simpa.utils.calculate</code> (module)	<code>STRUCTURE_FIRST_EDGE_MM</code> (simpa.utils.tags.Tags attribute)
<code>simpa.utils.constants</code> (module)	<code>STRUCTURE_RADIUS_MM</code> (simpa.utils.tags.Tags attribute)
<code>simpa.utils.deformation_manager</code> (module)	<code>STRUCTURE_RADIUS_VARIATION_FACTOR</code> (simpa.utils.tags.Tags attribute)
<code>simpa.utils.dict_path_manager</code> (module)	<code>STRUCTURE_SECOND_EDGE_MM</code> (simpa.utils.tags.Tags attribute)
<code>simpa.utils.libraries</code> (module)	<code>STRUCTURE_SEGMENTATION_TYPE</code> (simpa.utils.tags.Tags attribute)
<code>simpa.utils.libraries.literature_values</code> (module)	<code>STRUCTURE_START_MM</code> (simpa.utils.tags.Tags attribute)
<code>simpa.utils.libraries.molecule_library</code> (module)	<code>STRUCTURE_THIRD_EDGE_MM</code> (simpa.utils.tags.Tags attribute)
<code>simpa.utils.libraries.spectra_library</code> (module)	<code>STRUCTURE_TYPE</code> (simpa.utils.tags.Tags attribute)
<code>simpa.utils.libraries.structure_library</code> (module)	<code>STRUCTURE_X_EXTENT_MM</code> (simpa.utils.tags.Tags attribute)
<code>simpa.utils.libraries.tissue_library</code> (module)	<code>STRUCTURE_Y_EXTENT_MM</code> (simpa.utils.tags.Tags attribute)
<code>simpa.utils.path_manager</code> (module)	<code>STRUCTURE_Z_EXTENT_MM</code> (simpa.utils.tags.Tags attribute)
<code>simpa.utils.settings</code> (module)	<code>Structures</code> (class in <code>simpa.utils.libraries.structure_library</code> )
<code>simpa.utils.tags</code> (module)	<code>STRUCTURES</code> (simpa.utils.tags.Tags attribute)
<code>simpa.utils.tissue_properties</code> (module)	
<code>SIMPA_OUTPUT_NAME</code> (simpa.utils.tags.Tags attribute)	
<code>SIMPA_OUTPUT_PATH</code> (simpa.utils.tags.Tags attribute)	
<code>SIMPAJSONSerializer</code> (class in <code>simpa.io_handling.serialization</code> )	
<code>SIMPASerializer</code> (class in <code>simpa.io_handling.serialization</code> )	
<code>simulate()</code> (in module <code>simpa.core.simulation</code> )	
<code>SIMULATE_DEFORMED_LAYERS</code> (simpa.utils.tags.Tags attribute)	
<code>SIMULATION_PATH</code> (simpa.utils.tags.Tags attribute)	

subcutaneous\_fat()  
(simpa.utils.libraries.tissue\_library.TissueLibrary  
method)

## T

Tags (class in simpa.utils.tags)

TIME\_SERIES\_DATA (simpa.utils.tags.Tags attribute)

TIME\_SERIES\_DATA\_NOISE (simpa.utils.tags.Tags  
attribute)

TIME\_STEP (simpa.utils.tags.Tags attribute)

TISSUE\_PROPERTIES\_OUPUT\_NAME  
(simpa.utils.tags.Tags attribute)

TissueLibrary (class in  
simpa.utils.libraries.tissue\_library)

TissueProperties (class in  
simpa.utils.tissue\_properties)

to\_settings()  
(simpa.utils.libraries.structure\_library.Background  
method)

(simpa.utils.libraries.structure\_library.CircularTubularStructure  
method)

(simpa.utils.libraries.structure\_library.EllipticalTubularStructure  
method)

(simpa.utils.libraries.structure\_library.GeometricalStructure  
method)

(simpa.utils.libraries.structure\_library.HorizontalLayerStructure  
method)

(simpa.utils.libraries.structure\_library.ParallelepipedStructure  
method)

(simpa.utils.libraries.structure\_library.RectangularCuboidStructure  
method)

(simpa.utils.libraries.structure\_library.SphericalStructure method)

(simpa.utils.libraries.structure\_library.VesselStructure method)

TOTAL\_TIME (simpa.utils.tags.Tags attribute)

TUKEY\_WINDOW\_ALPHA (simpa.utils.tags.Tags  
attribute)

## U

UniformNoiseProcessingComponent (class in simpa.co  
re.processing\_components.noise.uniform\_noise)

UNITS\_ARBITRARY (simpa.utils.tags.Tags attribute)

UNITS\_PRESSURE (simpa.utils.tags.Tags attribute)

update\_internal\_properties() (simpa.utils.libraries.molec  
ule\_library.MolecularComposition method)

update\_settings\_for\_use\_of\_model\_based\_volume\_cre  
ator() (simpa.core.device\_digital\_twins.devices.pa\_dev  
ices.ithera\_msot\_acuity.MSOTAcuityEcho method)

UPSAMPLED\_DATA (simpa.utils.tags.Tags attribute)

US\_GEL (simpa.utils.tags.Tags attribute)

## V

VesselStructure (class in  
simpa.utils.libraries.structure\_library)

view\_absorption\_spectra() (in  
simpa.utils.libraries.spectra\_library) module

VOLUME\_CREATION\_MODEL\_SETTINGS  
(simpa.utils.tags.Tags attribute)

VOLUME\_CREATOR (simpa.utils.tags.Tags attribute)

VOLUME\_CREATOR\_SEGMENTATION\_BASED  
(simpa.utils.tags.Tags attribute)

VOLUME\_CREATOR\_VERSATILE  
(simpa.utils.tags.Tags attribute)

VOLUME\_NAME (simpa.utils.tags.Tags attribute)

VolumeCreationModelModelBasedAdapter (class in si  
mpa.core.volume\_creation\_module.volume\_creation\_m  
odule\_model\_based\_adapter)

VolumeCreationModuleSegmentationBasedAdapter  
(class in simpa.core.volume\_creation\_module.volume\_  
creation\_module\_segmentation\_based\_adapter)

VolumeCreatorModuleBase (class in  
simpa.core.volume\_creation\_module)

## W

warning() (simpa.log.file\_logger.Logger method)

WAVELENGTH (simpa.utils.tags.Tags attribute)

WAVELENGTHS (simpa.utils.tags.Tags attribute)





# Python Module Index

## S

### simpa

[simpa.core.acoustic\\_forward\\_module](#)

[simpa.core.acoustic\\_forward\\_module.acoustic\\_forward\\_module\\_k\\_wave\\_adapter](#)

[simpa.core.device\\_digital\\_twins.devices.detection\\_geometries.curved\\_array](#)

[simpa.core.device\\_digital\\_twins.devices.detection\\_geometries.detection\\_geometry\\_base](#)

[simpa.core.device\\_digital\\_twins.devices.detection\\_geometries.linear\\_array](#)

[simpa.core.device\\_digital\\_twins.devices.detection\\_geometries.planar\\_array](#)

[simpa.core.device\\_digital\\_twins.devices.illumination\\_geometries.illumination\\_geometry\\_base](#)

[simpa.core.device\\_digital\\_twins.devices.illumination\\_geometries.pencil\\_array\\_illumination](#)

[simpa.core.device\\_digital\\_twins.devices.illumination\\_geometries.pencil\\_beam\\_illumination](#)

[simpa.core.device\\_digital\\_twins.devices.illumination\\_geometries.slit\\_illumination](#)

[simpa.core.device\\_digital\\_twins.devices.pa\\_devices.ithera\\_msot\\_acuity](#)

[simpa.core.device\\_digital\\_twins.devices.pa\\_devices.ithera\\_msot\\_invision](#)

[simpa.core.device\\_digital\\_twins.devices.pa\\_devices.ithera\\_rsom](#)

[simpa.core.optical\\_simulation\\_module](#)

[simpa.core.optical\\_simulation\\_module.optical\\_forward\\_model\\_mcx\\_adapter](#)

[simpa.core.processing\\_components](#)

[simpa.core.processing\\_components.noise.gamma\\_noise](#)

[simpa.core.processing\\_components.noise.gaussian\\_noise](#)

[simpa.core.processing\\_components.noise.poisson\\_noise](#)

[simpa.core.processing\\_components.noise.salt\\_and\\_pepper\\_noise](#)

[simpa.core.processing\\_components.noise.uniform\\_noise](#)

[simpa.core.reconstruction\\_module](#)

[simpa.core.reconstruction\\_module.reconstruction\\_module\\_delay\\_and\\_sum\\_adapter](#)

[simpa.core.reconstruction\\_module.reconstruction\\_module\\_delay\\_multiply\\_and\\_sum\\_adapter](#)

[simpa.core.reconstruction\\_module.reconstruction\\_module\\_signed\\_delay\\_multiply\\_and\\_sum\\_adapter](#)

[simpa.core.reconstruction\\_module.reconstruction\\_module\\_time\\_reversal\\_adapter](#)

[simpa.core.reconstruction\\_module.reconstruction\\_utils](#)

[simpa.core.simulation](#)

[simpa.core.volume\\_creation\\_module](#)

[simpa.core.volume\\_creation\\_module.volume\\_creation\\_module\\_model\\_based\\_adapter](#)

[simpa.core.volume\\_creation\\_module.volume\\_creation\\_module\\_segmentation\\_based\\_adapter](#)

[simpa.io\\_handling](#)

[simpa.io\\_handling.io\\_hdf5](#)

[simpa.io\\_handling.serialization](#)

[simpa.log](#)

[simpa.log.file\\_logger](#)

[simpa.utils.calculate](#)

[simpa.utils.constants](#)

[simpa.utils.deformation\\_manager](#)

[simpa.utils.dict\\_path\\_manager](#)

[simpa.utils.libraries](#)

[simpa.utils.libraries.literature\\_values](#)

[simpa.utils.libraries.molecule\\_library](#)

[simpa.utils.libraries.spectra\\_library](#)

[simpa.utils.libraries.structure\\_library](#)

[simpa.utils.libraries.tissue\\_library](#)

[simpa.utils.path\\_manager](#)

[simpa.utils.settings](#)

[simpa.utils.tags](#)

[simpa.utils.tissue\\_properties](#)